


CHM[®]

STERNO-COSTAL PLATE

- *IMPLANTS*
- *INSTRUMENT SET 40.5841.000*
- *SURGICAL TECHNIQUE*



SYMBOLS DESCRIPTION



Caution - pay attention to a special procedure.



Perform the activity under X-Ray control.



Information about the next stages of a procedure.



Proceed to the next stage.



Return to the specified stage and repeat the activity.

www.chm.eu

Document No ST/42A
Date of issue 01.02.2012
Review date P-005-23.06.2023

*The manufacturer reserves the right to introduce design changes.
Updated INSTRUCTIONS FOR USE are available at the following website: ifu.chm.eu*

1. INTRODUCTION	5
1.1. BACKGROUND INFORMATION	5
1.2. INDICATIONS	5
1.3. CONTRAINDICATIONS	5
1.4. WARNINGS AND PRECAUTIONS	5
2. IMPLANTS	6
3. INSTRUMENT SET	8
4. SURGICAL TECHNIQUE	9
4.1. PATIENT POSITION	9
4.2. SURGICAL APPROACH	9
4.3. IMPLANT'S SELECTION	9
4.4. PLATE SHAPING	10
4.5. TUNNEL PERFORMANCE FOR SEATING OF THE PLATE	10
4.6. INITIAL CORRECTION OF DEFORMATION	11
4.7. IMPLANT'S INSERTION	11
4.8. INITIAL EVALUATION OF CORRECTION	12
4.9. CROSSWISE PLATE ATTACHMENT	13
4.10. IMPLANT REMOVAL	13

1. INTRODUCTION

1.1. BACKGROUND INFORMATION

Sterno-Costal Plate is used in the treatment of the funnel chest deformity, so-called pectus excavatum, with use of dr Nuss method (*minimally invasive repair of pectus excavatum – MIRPE*). The method involves insertion of the plate under sternum in order to achieve deformity correction. The procedure is minimally invasive and is associated with reduced operating time and minimal blood loss.

The method ensures:

- excellent cosmetic results,
- proper growth of lungs and heart,
- improvement of chest elasticity,
- remarkable breathing improvement.

The average time needed for a patient to return to normal daily activities is about a month.

1.2. INDICATIONS

Sterno-costal plate is used for treatment of chest deformities, especially the funnel chest (*lat. Pectus Excavatum*). The method of treatment using sterno-costal plate is designed especially for growing children (*when the ribs and costal cartilage are neither too malleable nor too rigid*). Optimal recommended age for reconstruction ranges from 12 to 16 years.

1.3. CONTRAINDICATIONS

Do not use the sterno-costal plate for patients with:

- mental illness or neurological disease,
- insufficient bone and fibrous tissue strength,
- infection.

The above list is not exhaustive.

For further information on:



- **adverse effects,**
- **warnings,**
- **sterilization,**
- **pre- and post-operative recommendations,**

please refer to the Instructions for Use (IFU) for sterno-costal plate.

1.4. WARNINGS AND PRECAUTIONS

The surgeon should avoid bending in sharp curves, reverse bending, and bending the implant at a hole. As a result of inappropriate shaping, size selection, wrong stabilization and fixation and patient's non-compliance with recommendations regarding principles of behavior during postoperative period, displacement or rotation of the implant may occur. It may lead in consequence to damage of tissue or organs adjacent to the implant.

During the implantation procedure, extreme caution should be taken to avoid contact of the implant or instruments with heart and lungs, because it may lead in consequences to permanent damage of these organs or in extreme case - death of the patient. After achieving the stable correction of deformity, implant has to be removed. After implant removal patient should be followed by postoperative monitoring to check for reoccurrence of the deformity.




This surgical technique is intended as a guide only. The selection of surgical technique adequate for specific patient remains surgeon's responsibility.
The implantation shall be carried out by the surgeon familiar with adequate rules and operating techniques, who acquired practical skills of using ChM instrument set for ChM sterno-costal plate.

2. IMPLANTS


Sterno-costal plate

L [mm]	Tytan Titanium	
	B=2.5[mm]	B=3.0[mm]
180	3.6116.180	3.6124.180
205	3.6116.205	3.6124.205
230	3.6116.230	3.6124.230
255	3.6116.255	3.6124.255
280	3.6116.280	3.6124.280
305	3.6116.305	3.6124.305
330	3.6116.330	3.6124.330
355	3.6116.355	3.6124.355
380	3.6116.380	3.6124.380
405	3.6116.405	3.6124.405
430	3.6116.430	3.6124.430




Crosswise plate

L [mm]	Tytan Titanium
45	3.6118.045
50	3.6118.050
55	3.6118.055



Dual crosswise plate


L [mm]	A [mm]	Tytan Titanium
60	15	3.6119.015
65	20	3.6119.020
70	25	3.6119.025
75	30	3.6119.030
80	35	3.6119.035



⚠ The Crosswise plates [3.6118.xxx], [3.6119.xxx] are intended for use with Plate blocker [3.6117.000] only.

Plate-blocker

Tytan Titanium
3.6117.000




Palette for implants - Sterno-costal plates 40.5843.000




Crosswise plate 0°

Tytan Titanium
3.6114.000




Crosswise plate

	Tytan Titanium	
α	left	right
10°	3.6113.010	3.6112.010
20°	3.6113.020	3.6112.020



Dual crosswise plate


L [mm]	A [mm]	Tytan Titanium
60	15	3.6120.015
65	20	3.6120.020
70	25	3.6120.025
75	30	3.6120.030
80	35	3.6120.035



The Crosswise plates [3.6112.xxx], [3.6113.xxx], [3.6114.xxx], [3.6120.0xx] are intended for use with Plate-blocker [3.6115.000] only.

Plate-blocker









Tytan Titanium
3.6115.000



Palette for implants - Sterno-costal plates 40.5843.200



3. INSTRUMENT SET

40.5841.000	Name	Pcs	Catalogue no.
	Raspatory L=510mm	1	40.6142.510
	Bender for plates	1	40.5848.000
	Persuader	1	40.6143.000
	Persuader	1	40.6144.000
	Torque handle 2Nm	1	40.6657.000
	Screwdriver tip T15	1	40.5677.000
	Plate trial L=180mm	1	40.5844.180
	Plate trial L=205mm	1	40.5844.205
	Plate trial L=230mm	1	40.5844.230
	Plate trial L=255mm	1	40.5844.255
	Plate trial L=280mm	1	40.5844.280
	Plate trial L=305mm	1	40.5844.305
	Plate trial L=330mm	1	40.5844.330
	Plate trial L=355mm	1	40.5844.355
	Plate trial L=380mm	1	40.5844.380
	Plate trial L=405mm	1	40.5844.405
	Plate trial L=430v	1	40.5844.430
	Container lid 9x4	1	14.0911.102
	Tray for implants - Sterno-costal plates	1	14.0911.601
	Container 9x4H	1	14.0911.101

4. SURGICAL TECHNIQUE


4.1. PATIENT POSITIONING

The patient is placed on his back with both arms abducted at the shoulders (up to angle of 90°), and the forearms bent to a right angle which allows easy access to the lateral chest wall. Correct positioning of the arms prevent neurological complications.



4.2. SURGICAL APPROACH

Mark the deepest portion of the chest using a sterile marker (if the deepest portion of the chest is below the sternum, mark the deepest point on the sternum). Determine the intercostal space (on the both side of the chest) located in line with the point set on sternum (or passing close to that point). The entry point (transverse lateral incision) is performed on the extension of the drawn line, between anterior axillary line and midaxillary line.



4.3. IMPLANT'S SELECTION

Because of the numerous types of chest deformation, the selection of the implant (length) should be preceded by proper measurements. The shape and length of the plate is determined by shape of deformation. The appropriate length selection allows to gain suitable plate stability. The measurements has to be taken before the surgery and confirmed during the procedure. Use plate trial [40.5844.xxx].



NOTE:

The length of required Sterno-costal plate must be smaller than the measured distance because the trial measures the external dimension of the chest and the implant traverses interiorly



4.4. PLATE SHAPING

Plate shaping has to be made with use of the bender [40.5848.000], suitable for the shape of chest deformity and implant insertion site.



For sliding crosswise plates [3.6118] or [3.6119] (blue color), during shaping, their chamfered ends have to be directed outside (IV.4a).

However, for crosswise plates applied on top surface of the sterno-costal plate [3.6112], [3.6113] or [3.6114] (gold color coded), during shaping, the chamfered ends of the sterno-costal plate have to be directed inside (IV.4b).

NOTE: As a result of acting forces deriving from the sternum and from the pressure present in the chest, the increase of the curvature in the middle of plate might be necessary in order to level the initial deformation of the plate caused by above mentioned forces.



Contouring of Sterno-Costal Plate should only be done with proper equipment.

The operating surgeon should avoid:



- reverse bending of the implant when contouring; reverse bending causes surface defects and internal stresses which may significantly decrease the fatigue life and may result in potential fracture of the implant;
- sharp bending of the implant (bending of a short segment and/or with a small bending radius);
- bending of the implant at a hole.

Do not excessively bend the ends of the plate in portions of the crosswise plates seating and the locking holes. Excessive bending may lead to deformation of the locking thread and/or may cause difficulties with proper fixation of the crosswise plates.



4.5. TUNNEL PERFORMANCE FOR IMPLANTATION OF THE PLATE




During the chest surgery, extreme caution should be taken. Contact of the implant or instruments with heart and lungs may lead in consequences to permanent damage of these organs or in extreme case - death of the patient. Therefore, to increase safety and facilitate the procedure of plate insertion it is needed to use the thoracoscope to visualize the chest organs.

Perform 2.5cm long incisions at both sides of the chest along indicated incision lines.




Choose the preparer size to fit the patient's chest.

Insert the raspatory [40.6142.510] to the made incision from the right side of the chest. With gentle movements, push the raspatory across mediastinum just below the sternum (*the tip of the preparer should stay in continuous contact with the sternum*), making the tunnel for the implant placement. At the final stage, the tip of the raspatory has to be advanced through the opposite incision in the intercostal space.




4.6. INITIAL CORRECTION OF DEFORMATION

Make the initial deformity correction by lifting both ends of the raspatory and pressing the places above and below the sternum in order to stretch the connective tissues. The initial correction facilitates the subsequent rotation of the plate and increases the stability of implant fixation.




4.7. IMPLANT INSERTION


Fix the umbilical tape to the end of the raspatory, and then withdraw the instrument, dragging the end of the tape to the other side of the chest.



Fix the plate to the end of the tape (at the right side of the chest). Gently pull the plate through the tunnel made earlier (the convexity of the plate should be faced down).



Once the plate is positioned, use persuader [40.6143.000] to turn the plate by 180° (directing the curved ends downwards), causing elevation of the sternum and correction of deformity. In order to facilitate the plate rotation, the second persuader may be used [40.6144.000].



4.8. INITIAL EVALUATION OF CORRECTION

Initial evaluation of correction is aimed to define what kind of plate stabilization will be required for its stable fixation – the usage of one or two crosswise plates decreases the probability of plate rotation.





It is recommended to use one crosswise plate. The usage of the second crosswise plate can be determined by patient's age, physical activity, muscle building and it depends on surgeon's decision.

For insufficient correction (elderly patients, patients with deeper deformity), it may be necessary to place the second sterno-costal plate. Additional plate is placed above or below the first plate.

4.9. CROSSWISE PLATE INSTALLATION

Fit the crosswise plate on the chosen end of the sterno-costal plate and then determine its position (*the most common is 1-2cm from the end of the sterno-costal plate*). In order to fix the gold color coded crosswise plate, the plate blocker [3.6115.000] shall be used, whereas to fix the blue color coded crosswise plates, one shall use the plate blocker [3.6117.000].



It is also possible (*optionally*) to fix position of the sliding crosswise plate (*blue color coded*) on the sterno-costal plate by means of multiple strong figure-8 sutures.

Then, by means of strong sutures, fix the implants to the chest wall muscles, using the holes in the sterno-costal plate and crosswise plate.

Before the wound closure, place the patient in Trendelenburg position, inflate the lungs with air and give positive end-expiratory pressure (*PEEP*) to prevent pleural air trapping. Cover the implants with surrounding soft tissue and skin and then close the wound with an absorbing suture and a dressing.



A chest radiograph should be obtained postoperatively due to possibility of pneumothorax occurrence and to confirm the proper position of the implant.

4.10. IMPLANT REMOVAL

The procedure of the plate removal is performed with total anesthesia in ambulatory conditions. During the procedure, patient is placed in supine position with arms abducted. In order to remove the plate, incisions are made in the same locations as during the implantation, allowing access to the plate, crosswise plates and sutures. The plate is withdrawn (*after earlier suture removal*) pulling the one end through the incision and turning the patient in the opposite direction. After the implant removal, the wound is closed by means of absorbable sutures. Postoperative chest radiograph is recommended.

ChM sp. z o.o.

Lewickie 3b
16-061 Juchnowiec Kościelny
Poland

tel. +48 85 86 86 100

fax +48 85 86 86 101

chm@chm.eu

www.chm.eu



CE 0197