CHM[®]

30.0008÷30.0019

PL	OPASKA ZACISKOWA	3
EN	TOURNIQUET	10
RU	ТУРНИКЕТ	17

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Tab. 1. Description of available tourniquets

9. INDICATIONS FOR USE

Tourniquet, in combination with an external source of compressed air or nitrogen, is used to control blood flow in the limb and produce a bloodless surgical field. Tourniquet control unit or hand pump (*compressed air*) produced by **ChM** may be examples of external source of compressed air or nitrogen.

10. CONTRAINDICATIONS

- Open fracture of the lower limb.
- Post-traumatic hand reconstruction.
- Severe crush injuries.
- Elbow surgery (*with excessive swelling*).
- Severe hypertension.
- Skin grafts in which all bleeding points must be distinguished.
- Compromised circulation e.g. peripheral arterial disease.
- Diabetes.
- Sickle cell disease is a relative contraindication against use of a tourniquet (*See PRECAUTIONS*).
- Avoid using a tourniquet in patients undergoing secondary or delayed treatment after fixation.

NOTE:

The abovementioned list of contraindications is not exhaustive.

11. USE

Depending on the area of application and the type of surgical procedure, tourniquets of different structure may be used so that they are more effective. **ChM** offers several types of tourniquets differing in shape, length and width, allowing the surgeon to choose a tourniquet appropriate to the location and type of procedure to be performed.

11.1. PLACE OF USE

The tourniquet shall be applied on a limb at the place of a large muscle mass which helps to protect vessels and nerves from damage.

The typical location of the tourniquet:

- arms and thighs - where the limb has the the largest circumference closest to the operation field
- forearm - the central part of the forearm
- shank - the proximal edge of the tourniquet should be placed on the largest circumference of the calf
- foot - 1/3 of the lower part of the shank with the distal edge of the tourniquet passing above the ankle

11.2. SINGLE AND DUAL TOURNIQUETS

Depending on the surgical procedure, single or dual tourniquets may be used. Single tourniquets are mostly used with general anesthesia, brachial plexus blocks and central blocks.

Dual tourniquets are mostly used with local intravenous anesthesia (*also called Bier's block*). Thanks to the independent inflation of the proximal and distal section of the tourniquet, the surgical procedure is more patient-friendly, especially during longer operations.

11.3. LENGTH

The length of a tourniquet is chosen on the basis of the length of the limb circumference measured at the point where the tourniquet is to be applied. Knowing the circumference of the limb and the place of application, select the tourniquet from the table (Tab. 1).

**NOTE:**

Use of excessively short or long tourniquets can lead to its unstable positioning on a limb, problems with full occlusion of blood vessels and damage to the skin.

11.4. WIDTH

As wide as possible tourniquets should be applied. A wider tourniquet requires lower occlusion pressure than a narrow one what reduces the risk of complications associated with its use.

11.5. TOURNIQUET PRESSURE

Tourniquet pressure used during surgery depends on many variables which include: the patient's age, skin condition, a place the tourniquet is applied, the shape and size of the tourniquet. According to the publicly available literature data, occlusion pressure:

- for the upper limb should be about 75÷100mmHg above systolic blood pressure,
- for the lower limb should double the value of systolic blood pressure.

11.6. PREPARING A PATIENT FOR SURGERY

Each patient that is scheduled for surgery with the use of a tourniquet, should undergo a preoperative evaluation including, inter alia: determination of the physical condition of the patient, a medical history and medication review. On the basis of the preoperative evaluation, possible complications that may occur during the surgery are determined and whether the patient is allowed to be treated with the use of the tourniquet. The measurement of the blood pressure and the length of the limb at the point where the tourniquet will be applied should also be performed preoperatively.

11.7. TOURNIQUET USE

Carefully check the condition and tightness of the tourniquet prior to each use. Should a defect or doubts as to the condition of the tourniquet appear, it is unacceptable to continue the use of the tourniquet. Service and warranty repairs may be performed only by the manufacturer - **ChM**.

Apply disposable protective cuff on the selected preoperatively tourniquet that would protect the tourniquet from blood, body fluids or other infectious agents contamination. The skin under the tourniquet should be protected from direct exposure to that tourniquet by applying padding underneath, keeping in mind though that excessively thick padding may increase occlusion pressure. Moreover, do not allow any fluids to get under the tourniquet because they can cause skin injuries. When applying the tourniquet, make sure that the surfaces that come into direct contact with the padding (*skin*) were not wrinkled because it can lead to skin damage after inflating the tourniquet. Applied tourniquet should adhere to the limb with its entire perimeter. The tourniquet feeder should be placed outside of the limb and operating field. Avoid accidental pressure of the feeder on nerves or its bend. Connect the tourniquet to a tourniquet control unit and start inflating until the established occlusion pressure is reached. Inflating the tourniquet should proceed quickly and smoothly so that arteries and veins

are closed at the same time.

As a general accepted rule, the maximum time of the tourniquet being inflated (*without its deflating*), in the case of a healthy adult patient, should not exceed 120 minutes. During the surgery, to reduce the risk of complications, intraoperative control of blood pressure, the pressure in the tourniquet and the time the tourniquet is applied should be conducted.

Before deflating the tourniquet, it is recommended to close the wound. Deflating the tourniquet before wound closure increases blood loss. Deflating the tourniquet should proceed quickly and smoothly. After deflating, immediately remove the tourniquet and the padding from the limb so that they will not block the return of venous blood and resultant blood accumulation in the operating field.

12. SIDE EFFECTS

After tourniquet use, dull ache in the limb may appear.

Prolonged ischemia caused by use of the tourniquet can lead to temporary or permanent damage to tissues, blood vessels and nerves and can cause changes in blood clotting.

Excessively high pressure in the tourniquet can cause paralysis of limbs.

Excessively low pressure in the tourniquet can cause passive congestion of limbs with possible irreversible loss of function.

Intraoperative bleeding may be caused by:

- blood remaining in the limb due to insufficient emptying the limb of blood
- the use of improper occlusion pressure
- blood penetrating from the nourishing vessels of the long bones.

13. CONSTRUCTION

The tourniquet is made of polyamide fabric.

The latex bladder connected with feeder is placed in the tourniquet pocket.

The quick-fit port connecting the tourniquet to the tourniquet control unit is placed at the end of the feeder. The tourniquet is reinforced by silicone insert to obtain greater stiffness.

There are special Velcro and protective tapes on the inner side of the cuff to protect the tourniquet against loosening when applied on the limb and inflated.

14. DISINFECTION AND CLEANING

The tourniquet should be cleaned by hand with a soft brush dampened in warm water and a small amount of soap causing no foaming, unless the tourniquet got contaminated with infectious agents, e.g. blood. Only exterior surfaces of the rubber feeder and cuff should be cleaned, avoiding the excessive soaking of the tourniquet material. It is not allowed to: immerse the tourniquet in the water or aqueous solutions of washing and disinfecting agents, machine-wash the tourniquet or to pull out the inner latex bladder from the cuff. Wipe the washed tourniquet with a damp lint-free cloth and air-dry (*do not roll*). When dried, the tourniquet may be wiped with cloths soaked in isopropyl alcohol.

Should the tourniquet have contact with infectious agents, e.g. blood, for cleaning and disinfection use mild washing and disinfecting agents approved for use with medical devices. Follow the instructions provided by the manufacturer of these agents. It is recommended to use aqueous solutions of washing or washing-disinfecting agents with a pH value between 7 and 10.8.

1. Remove blood and/or dirt from the product using disposable cloths or paper towels.
2. Prepare a fresh detergent solution using e.g., Neodisher MediClean Forte from Dr.Weigert at a temperature of $40^{\circ}\text{C} \pm 2$, using 15ml of the agent per 1l of water. Thoroughly clean the surface of the cuff with a soft, lint-free cloth.
3. Dry the device thoroughly using disposable, soft, lint-free cloth.
4. Prepare an aqueous solution of disinfectant, e.g. Neodisher Septo Active from Dr.Weigert of 2% concentration (*20g of agent per 1 liter of water*) and temperature $20^{\circ}\text{C} \pm 2$. Thoroughly clean the surface of the cuff with a soft, lint-free cloth.
5. After disinfecting, clean the tourniquet with a damp cloth to remove residual solution that may have an allergic and irritating effect and shorten the life of the product. Dry the tourniquet in the open air - do not roll.
6. Clean the tourniquet with cloths soaked in isopropyl alcohol.

NOTE:

The tourniquet, when used in surgery, should always be placed in a protective cuff which protects it from contamination by infectious agents.

15. STORAGE

The tourniquet shall be stored unwound or slightly twisted in a roll, in a dry, not in direct sunlight place. The tourniquet storage should be preceded by its washing and thorough drying.

16. WARRANTY AND SERVICE

The manufacturer provides annual warranty for the usage of the tourniquet provided that the above principles are observed.

Unauthorised changes introduced by the user might render the warranty void.

ChM sp. z o.o. does not specify the maximum number of cycles for use of tourniquets. The tourniquet life is dependent on many factors, including the duration of each use, the pressures used, re-processing.

Each time before re-use, inspect the functionality of the product. In case of damage, the tourniquet must not be used for surgery.

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