ChM®

ChM Power Tools
CharDrive *system*

PL	SYSTEM NAPIĘDU MEDYCZNEGO CharDrive	3
EN	CharDrive MEDICAL DRIVE SYSTEM	39
RU	СИСТЕМА МЕДИЦИНСКОГО ПРИВОДА CharDrive	75

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1. INTRODUCTION	40
1.1. General information	40
2. OPERATING THE CharDrive SYSTEM	44
2.1. Description of CharDrive elements	44
2.2. CharDrive Medical drive system in use	48
3. CLEANING, DISINFECTION, STERILIZATION	57
3.1. General information	57
3.2. Handpiece and attachments	58
3.3. Storage and handling before use	61
3.4. Batteries and battery charger	62
4. REPAIRS	63
5. LIFETIME	63
6. DISPOSAL	63
7. TROUBLESHOOTING	64
7.1. CharDrive Handpiece	64
7.2. CharDrive Battery	65
7.3. Attachments and cutting elements	65
7.4. CharDrive battery charger	66
8. TECHNICAL DATA	67
8.1. Duty cycle	67
8.2. CharDrive Handpiece	68
8.3. CharDrive Battery	68
8.4. CharDrive Battery charger	68
9. ELECTROMAGNETIC COMPATIBILITY IN ACCORDANCE WITH EN 60601-1-2	69
9.1. CharDrive Handpiece	69
9.2. CharDrive Battery charger	71
10. EXPLANATION OF SYMBOLS USED	73

Prior to first use, read the instructions for use.

1. INTRODUCTION

1.1. GENERAL INFORMATION

1.1.1. INTENDED USE

The **CharDrive** Handpiece is a battery-powered medical device (*powered by a lithium-ion battery*) intended for use in general surgical procedures including drilling and reaming of bone tissues. The **CharDrive** Medical Drive System consists of: handpiece, 2 rechargeable batteries, battery charger, cover and a set of attachments.

1.1.2. SAFETY INSTRUCTIONS⁽¹⁾

The **CharDrive** Medical Drive System can only be used after careful consultation of the system's Instructions for Use. It is recommended that an alternative (*spare*) set of devices be available during the surgery since unexpected technical problems can occur.

The **CharDrive** Medical Drive System can only be used by physicians and trained medical personnel. It is recommended to take utmost care and to avoid haste when using the system. Check all components of the system prior to use. Lack of basic knowledge about the functionality of components can result in injury to a patient, medical staff or operator.

Verify the system components and cutting elements for their physical properties before use. DO NOT use any of the components if damaged.

DO NOT use this system in presence of oxygen, nitrous oxide or a mixture consisting of flammable anesthetic and air.

To ensure operation of the system, use only **ChM** original accessories. Using accessories of another manufacturers may result in unplanned functioning of the device and loss of warranty rights. DO NOT modify any system components yourself. Failure to follow the above-mentioned can result in injury to a patient, medical staff or operator.

Before first and every use, the **CharDrive** Medical Drive System (*except for battery charger and batteries*) has to be re-processed thoroughly. Protective covers and foils must be fully removed before sterilization.

For the **CharDrive** Medical Drive System to function properly, **ChM** recommends that the system be cleaned after each use in accordance with the instructions provided in the chapter CLEANING, DISINFECTION, STERILIZATION.

Following the recommendations provided can considerably extend the serviceable life of the system.

Successful surgery depends on the condition of the instruments used. Therefore, verify the condition of the used cutting elements after every use for wear and/or damage and replace them if necessary. It is recommended

⁽¹⁾ Should the above-mentioned information be insufficient, please, contact the manufacturer.

to use new **ChM** cutting instruments for every surgery. Cutting instruments, while in use, must be cooled with irrigation to prevent tissues heat necrosis.

Vibrations of the attachments, components cooperating with the handpiece, can cause, e.g., tissue or bone injuries, oversized drilled holes for locking screws, excessive reaming of the medullary canal, or incorrect insertion of e.g. Kirschner wire. Should vibration occur, refer to the TROUBLESHOOTING section.

DO NOT use excessive pressure when working with cutting elements to prevent their breakage/ damage. Excessive pressure, especially during high rotation, can cause significant bending of the cutting elements which can consequently lead to: tissue damage, handpiece being torn from the operator's hand or cutting instrument's elements being spread over long distance. Failure to follow the above-mentioned can result in injury to a patient, medical staff or operator.

DO NOT reuse single-use cutting elements. Failure to follow the above-mentioned can result in injury to a patient, medical staff or operator.

DO NOT use a product that has been in contact with tissues and body fluids of another patient without prior washing, disinfection and sterilization of the product, since the potential risk for cross infection that may include viruses, bacteria, and prions occurs.

The user of the product is fully responsible for the proper use of the instruments during the surgery.

When not in use, before removing or installing co-operating accessories, and when handed over to another user, the handpiece should always be switched into safety mode. Failure to follow the above-mentioned can result in injury to a patient, medical staff or operator.

When not in use, the handpiece should be put on its side to prevent it from falling to the ground, should instability occur.

Before using the handpiece, make sure that the attachment is properly secured to the handpiece and that the cutting elements are properly installed in the attachment, e.g. by operating the set in the air. Failure to follow the above-mentioned can result in injury to a patient, medical staff or operator.

DO NOT put your hand on or hold any rotating part while using the handpiece. Failure to follow the above-mentioned can result in injury to a patient, medical staff or operator.

DO NOT overload the handpiece. Failure to follow the above-mentioned can result in motor and/ or battery damage.

If the handpiece is used in implantation procedures in combination with implant systems, refer to the appropriate *"Operating Technique"* for this implant system.

DO NOT dismantle or service the product yourself. The product should be sent to the manufacturer for repair. Failure to follow the above-mentioned can result in electric shock or fire.

DO NOT plug in or use the battery charger with a voltage other than that indicated on the back panel of the unit.

For batteries charging, use always original battery charger. Failure to follow the above-mentioned can result in injury to a patient, medical staff or operator.

DO NOT use any battery charger accessories that are not recommended or sold by **ChM**. Failure to follow the above-mentioned can result in injury to a patient, medical staff or operator.

DO NOT use the battery charger with a damaged cord or plug.

The power cord of the battery charger should be connected to the power supply socket equipped with protective grounding.

Prior to cleaning, ALWAYS unplug the battery charger from the power supply to minimize the risk of electric shock.

For information on electromagnetic compatibility (EMC), go to section ELECTROMAGNETIC COMPATIBILITY of this Instructions for Use.

The application part of the product has a BF-type of protection against electric shock. This product is suitable for use in patients according to EN 60601-1.

To ensure proper functioning of the handpiece, **ChM** recommends its annual inspection (*applies to a device for which the expected lifetime has not yet been exceeded*). Only trained and experienced personnel of the manufacturer can perform such inspections. For torque accessories, the calibration date must be strictly observed. The manufacturer shall not be liable for any damage resulting from unauthorized repairs of the product.

Before using the handpiece, make sure that the battery cover securing the battery in the handpiece is properly tightened.

The handpiece can only be operated with a fully charged battery. To do this, make sure that the battery is charged at the right time. It is recommended to recharge the battery used immediately after the treatment. Due to the fact that the batteries are consumables, it is recommended that they be replaced at least once a year.

The battery and battery charger must not be sterilized, washed, rinsed or soaked. These processes can destroy these devices and may cause secondary damage.

1.1.3. ACCESSORIES/SCOPE OF DELIVERY

The **CharDrive** Medical Drive System consists of:

- handpiece,
- 2 rechargeable batteries,
- battery charger,
- cover,
- a set of attachments to be installed in the handpiece.

For batteries charging use only the corresponding battery charger **[40.8045.500]**.

For the system to operate properly only **ChM** cutting instruments should be used.

The following components are essential to ensure proper functioning of the system:

- **CharDrive** Handpiece - Catalogue No. **[40.8000.400]**,
- **CharDrive** Battery - Catalogue No. **[40.8005.200]**,
- **CharDrive** Battery charger - Catalogue No. **[40.8045.500]**,
- **CharDrive** Cover - Catalogue No. **[40.8009.400]**,
- At least one attachment:
 - **CharDrive** AO quick coupling - Catalogue No. **[40.8010.200]**,
 - **CharDrive** Keyless attachment holder - Catalogue No. **[40.8014.200]**,
 - **CharDrive** Kirschner wire attachment - Catalogue No. **[40.8013.200]**,
 - **CharDrive** DIN reamer attachment - Catalogue No. **[40.8012.200]**,
 - **CharDrive** Hudson-type attachment - Catalogue No. **[40.8017.200]**,
 - **CharDrive** Stryker/Zimmer/Trinkle-type attachment - Catalogue No. **[40.8015.200]**,
 - **CharDrive** Harris-type attachment - Catalogue No. **[40.8016.200]**,
 - **CharDrive** AO-S-type attachment - Catalogue No. **[40.8018.200]**,
 - **CharDrive** Torque quick coupling 3/16 Square type-2Nm - Catalogue No. **[40.8022.220]**.

For sterilization of the elements of **CharDrive** Medical Drive System, use:

- **CharDrive**-sterilization stand - Cat. No. **[40.8041.000]**

The date of the next calibration, which is permanently marked on the torque instruments, must be strictly observed. Use of the torque instrument with an overdue calibration date may cause damage to the implant, instrument or loss of correction. Should some malfunction of the device be detected before the calibration date, e.g. as a result of intensive use, the instrument must immediately be sent to the manufacturer for calibration.

1.1.4. WARRANTY

The warranty for the device does not cover damage of any kind resulting from:

- improper use,
- improper in-hospital storage and transportation.


The manufacturer is not liable for damage resulting from repairs carried out by unauthorized persons.

2. OPERATING THE CharDrive SYSTEM

2.1. DESCRIPTION OF CharDrive ELEMENTS

2.1.1. CharDrive HANDPIECE [40.8000.400]


1. Locking sleeve for attachments - used to install an attachment to the handpiece; turn the sleeve clockwise.
2. Trigger for counterclockwise rotation - used to give the attachment counterclockwise rotation; press the trigger.
3. Trigger for clockwise rotation - used to give the attachment clockwise rotation; press the trigger.
4. Safety switch - depending on the switch position, the corresponding triggers may be locked or released.
 - a) the switch pressed from right - only clockwise rotation is available (*counterclockwise rotation is locked*)
 - b) the switch pressed from left - both triggers are locked - the safety mode is activated
 - c) neutral switch position - the clockwise and counterclockwise rotation is available; when two triggers are pressed simultaneously, the oscillation of an attachment is activated.
5. Battery cover - protects the battery from falling out from the handpiece - tighten the cover.
6. Information marked on the handpiece.
7. Mode change collar of the handpiece (*DRILL for drilling and REAM for reaming*) - turn the collar.
8. LED system that informs the user about the charge level of the battery at which the handpiece is currently operating.

A description of the symbols used can be found in chapter USED SYMBOLS of this Instructions for Use.

2.1.2. CharDrive BATTERY [40.8005.200]


1. Battery.
2. Handle for removing the battery from the handpiece.
3. Information marked on the battery.

A description of the symbols used can be found in chapter USED SYMBOLS of this Instructions for Use.

2.1.3. CharDrive BATTERY CHARGER [40.8045.500]

1. Battery charger ports with LED lights (for battery charging).
2. Graphic display.
3. Power button.
4. Power cord socket.
5. The vent.
6. Fuse socket.
7. Service socket - the socket is for service only, and it is forbidden to remove the plug and perform any operation during normal use of the charger.
8. Information marked on the back of the housing.
9. Information marked on the side of the housing.

A description of the symbols used can be found in chapter USED SYMBOLS of this Instructions for Use.

2.1.4. CharDrive COVER [40.8009.400]

The cover is used when the battery is placed in the handpiece under sterile conditions. The cover allows the non-sterile battery to be separated from the sterile outside of the handpiece casing.

1.1.1. ATTACHMENTS

	Attachment name	Cat. no.
	CharDrive AO quick coupling	40.8010.200
	CharDrive Keyless attachment holder	40.8014.200
	CharDrive DIN reamer attachment	40.8012.200
	CharDrive Hudson-type attachment	40.8017.200
	CharDrive Kirschner wire attachment	40.8013.200
	CharDrive Stryker/Zimmer/Trinkle-type attachment	40.8015.200
	CharDrive Harris-type attachment	40.8016.200
	CharDrive AO-S-type attachment	40.8018.200
	CharDrive Torque quick coupling 3/16 Square type-2Nm	40.8022.220

2.2. CharDrive MEDICAL DRIVE SYSTEM IN USE

2.2.1. INSERTION OF THE BATTERY INTO THE HANDPIECE

To ensure sterility, the battery should be inserted into the sterile housing of the handpiece by two people, one of whom is wearing sterile clothing:

1. The person wearing the sterile clothing holds the open, sterile handpiece with the open side up.
2. The person wearing the sterile clothing places the cover on the handpiece and makes sure it is correctly tightened. The cover ensures that the non-sterile battery does not have contact with the outside of the sterile handpiece.
3. The other person not wearing sterile clothing carefully inserts the non-sterile battery through the cover into the handpiece and presses firmly onto the battery to ensure that the battery is correctly seated in the handpiece. During insertion, it is important to ensure that the battery is properly positioned and that the person not wearing sterile clothing does not touch the outside of the sterile handpiece.
4. The person wearing non-sterile clothing grasps the cover and removes it from the handpiece.
5. The person wearing the sterile clothing places the battery cover on the handpiece and tightens it.

Fig. 1. Handpiece with the attached cover and inserted battery.

The handpiece can only be operated with a fully charged battery. To do this, make sure that the battery is charged at the right time. It is recommended to recharge the battery used immediately after the treatment. Due to the fact that the batteries are consumables, it is recommended that they be replaced at least once a year.

To ensure aseptic conditions, the battery must not be removed from the handpiece until the end of surgery. The battery has enough capacity to power the handpiece for the entire surgery. Additionally, the handpiece is equipped with a LED system that informs the user of the charge state of the battery on which the handpiece is currently operating.

Cover should be sterilized prior to first and after each use to ensure aseptic conditions when placing the non-sterile battery into the sterile handpiece.

2.2.2. REMOVAL OF THE BATTERY FROM THE HANDPIECE

To remove the battery from the handpiece, remove the battery cover, and holding the battery handle, remove the battery from the handpiece.

The handpiece must be held upright (Fig. 2) to prevent the battery from falling onto the floor.

Fig. 2. Position of the handpiece while inserting/ removing the battery.

2.2.3. CHARGING, STORING AND DISPOSAL OF THE BATTERY

Battery (*lithium-ion type*) is an element of the system that, like the battery charger, cannot be reprocessed (*sterilized*). To ensure the proper functioning of the product and to prolong its serviceable life, the following points must be observed:

2.2.3.1. Charging

- Use only the **ChM** battery charger to charge the batteries.
- The battery must be recharged immediately after use.
- Before charging the battery, make sure the charging socket is dry.
- Prior to use, the battery must be fully recharged even if it has not been used before.
- Charge the battery within an ambient temperature range from 10°C to 40°C.
- The charging cycle of a fully discharged battery takes up to 90 minutes.
- Depending on the procedure, the expected battery life is approximately 20 minutes (*minimum*).
- Use only batteries that have been charged in the last 48 hours.

2.2.3.2. Storage

- DO NOT expose the battery to temperatures over 50°C for time longer than 72h.
- DO NOT store the battery in the handpiece where it may be discharged.
- When stored off the drive, the battery will also gradually discharge and must be recharged at least once a year.

DO NOT:

- wash the battery,
- rinse the battery,
- sterilize the battery,
- drop the battery.

- Use only the ChM battery charger to charge the ChM batteries. The use of other chargers can damage the battery.
- Do not use faulty battery. The damaged battery should be sent to ChM.
- Use only the batteries in the designated handpiece.

2.2.3.3. Disposal

The handpiece is powered by the lithium-ion battery.

The batteries should be disposed of in an environmentally friendly manner. The batteries have been designed to be recycled as a whole without end user intervention. Disassembly is not required. Prior to recycling, contaminated products should be decontaminated. The European Battery Directive 2006/66/EC is applicable.

Symbol definition

Rechargeable lithium-ion battery.

In accordance with Directive 2006/66/EC, batteries must be collected and brought into use in accordance with the rules in force in the EU Member States.

2.2.4. OVERHEATING PROTECTION

Electrical systems can heat up if in constant use. To prevent overheating of the product, observe the handpiece duty cycles provided in the "Duty Cycle" section.

The thermal overload protection system prevents the handpiece from overheating.

If the handpiece becomes overheated during use, it automatically switches off and cannot be operated until "cooled down".

2.2.5. INSTALLING OF ATTACHMENTS AND CUTTING ELEMENTS TO THE HANDPIECE

Prior to assembly or disassembly of the attachment or mating accessories such as drill bits, the handpiece should always be set to the safety mode. Failure to follow the above-mentioned can result in injury to a patient, medical staff or operator.

Each time the attachment, cutting elements are attached, make sure that they have been properly secured to the handpiece by gently pulling on it and operating the set in the air.

Use only original ChM attachments and cutting elements.

Damage caused by the use of attachments and cutting elements of other companies are not covered by the manufacturer's warranty.

Use irrigation liquids to cool the cutting elements to prevent tissue necrosis.

Before and after surgery, check the condition of the cutting elements for wear and/ or damage and, if need be, replace them with new ones. It is up to the physician to decide on the use of the instrument. It is not allowed to use cutting elements with visible damage/ deformation. ChM recommends cutting instruments to be single use.

When using CharDrive Torque quick coupling 3/16 Square type-2Nm, the date of the next calibration, which is permanently marked on the instruments, must be strictly observed. Use of the torque instrument with an overdue calibration date may cause damage to the implant, instrument or loss of correction. Should some mal-function of the device be detected before the calibration date, e.g. as a result of intensive use, the instrument must immediately be sent to the manufacturer for calibration.

In order to install the selected attachment to the handpiece, turn the locking sleeve for attachments clockwise, first.

Then, holding the turned sleeve, place the attachment in the locking sleeve (Fig. 3) and press against the hand-

piece. Release the locking sleeve.
The attachment has been installed.

Attachments of ChM production have two shaped longitudinal protrusions for mounting in the handpiece. Adequately, the locking sleeve for attachments has two grooves that fit the protrusions and ensure correct mounting of the attachment (Fig.3).

After installing, make sure that the attachment has been properly secured to the handpiece by gently pulling on it and operating the set in the air.

Fig. 3. Attachment in the handpiece.

In order to change the attachment, stop the handpiece and switch into safety mode, first. Then, turn the locking sleeve for attachments clockwise and holding the turned sleeve, remove the attachment and place it gently and safely aside.

Install the selected attachment as described above.

CharDrive Keyless attachment holder

Keyless attachment holder is compatible with drills, reamers and Kirschner wires with a diameters from 0.3mm to 7.4mm.

To fit the cutting element in the keyless attachment holder, turn the front part of the attachment clockwise until an opening sufficient to accommodate the selected cutting element is obtained (Fig.4). Insert the selected cutting element and close the jaws of the attachment manually in a counterclockwise direction. After installing, make sure that the device has been properly secured to the attachment by gently pulling on it and operating the set in the air.

To remove the cutting element from the keyless attachment holder, turn the front part of the attachment clockwise until the instrument is loosened and remove the cutting element gently and safely.

ALWAYS manually close the jaws on the cutting element. Never close the jaws using the handpiece. Failure to follow the above-mentioned can result in injury to a patient, medical staff or operator.

CharDrive Kirschner wire attachment

To fit the Kirschner wire in the Kirschner wire attachment, place the selected wire (diameter between 0.8 and 3.2mm) in the attachment, and using the attachment lever, grasp the wire with adequate force. Holding the grip, make sure that the device has been properly secured to the attachment by gently pulling on it and operating the set in the air.

When inserting the wire into the bone, adjust the clamping force to the needs of the surgery (Fig.6).

To remove the cutting element from the attachment, release the lever and remove the wire gently and safely.

Fig.4. Open jaws of the keyless attachment holder

Fig.5. Closed jaws of the keyless attachment holder


Fig. 6. Kirschner wire attachment in use.

CharDrive: AO quick coupling; DIN reamer attachment; Hudson-type attachment; Stryker/Zimmer/Trinkle-type attachment; Harris-type attachment; AO-S-type attachment; CharDrive Torque quick coupling 3/16 Square type-2Nm

To fit the cutting element in the attachment, pull the collar of the attachment towards the handpiece, insert the compatible cutting element, and release the collar. The cutting element will automatically be mounted. After installing, make sure that the device has been properly secured to the attachment by gently pulling on it and operating the set in the air.

To remove the cutting element from the attachment, pull the collar towards the handpiece and remove the device gently and safely.

2.2.6. CharDrive BATTERY CHARGER IN USE

DO NOT plug in or use the battery charger with a voltage other than that indicated on the back panel of the unit.

For batteries charging, use always original battery charger. Failure to follow the above-mentioned can result in electric shock and/ or damage to the product.

DO NOT use any battery charger accessories that are not recommended or sold by ChM. Failure to follow the above-mentioned can result in electric shock and/ or damage to the product.

DO NOT use the battery charger with a damaged cord or plug.

DO NOT modify the grounding of the charger cable. The cable must be plugged in directly to the hospital power supply.

DO NOT dismantle or service the charger yourself. The product should be sent to the manufacturer for repair. Failure to follow the above-mentioned can result in electric shock or fire.

Always unplug the product from the power supply to minimize the risk of electric shock prior to cleaning processes of the charger.

For information on electromagnetic compatibility (EMC), go to section ELECTROMAGNETIC COMPATIBILITY of this Instructions for Use.

The battery and battery charger must not be sterilized, washed, rinsed or soaked. These processes can destroy these devices and may cause secondary damage.

2.2.6.1. Battery charger plugging in

WARNING: To avoid the risk of electric shock, the unit must only be connected to a power supply with protective grounding.

To connect the charger to a hospital power supply, connect the power cord to the socket on the back of the charger and then to the power supply socket with protective grounding.

Place the power cord of the charger so as to minimize the possibility of stepping on it, stumble over it or exposing it to other types of damage /stress.

DO NOT insert metallic elements into the charger, as the ports contacts may otherwise be damaged.

Before the first and every next use, the charger and its accessories (*such as power cord*) should be checked for possible damage to the casing, cable cracks, bent contacts, etc. If any defects are detected, the product cannot be used and should be returned to the manufacturer for repair purpose.

When charger is connected to the power supply, the start screen and batteries ports LED lighting (*blue*) appear (Fig. 1), and then 3 electronic batteries visualizations corresponding to 3 built-in charging ports will appear (Fig. 2).

Fig.1. Battery charger start screen

Fig.2. Information provided on the display

When inserting the battery into the charger port, its electronic equivalent displays the battery charge status (in %) and its charging process, illustrated by the increasing charging bars of the battery. Depending on the state of the battery discharge, the charging process will be visualized by four different colors displayed on the screen (Fig. 3).

Battery charge status (in%)	Color of battery charging process
0-30%	Red
35-60%	Orange
65-85%	Yellow
90-100%	Green

The charging cycle of a fully discharged battery takes up to 90 minutes. To ensure maximum operating time, prior to use, the battery must be recharged completely, even if it has not been used before. It is important to use only those batteries that have been charged in the last 48 hours.


Fig.3. Charging process and related LED port lighting

When the battery is 100% charged, the message 100% and 🟢 will appear on the display indicating that the battery is ready to operate. Green LED lighting will also flash (Fig. 4).

Occurrence of the message ❌ (Fig. 5) immediately after placing the battery in the charging port indicates a charging error. In this case, the charger and the batteries should be sent to the manufacturer to verify their functionality and possible repairs. A similar message will be displayed when the battery is removed from the charging port of the charger before the charging process is complete (*before the battery is fully charged*).

Fig.4. Information about correct, full charge of the battery and completion of charging process

Fig.5. Charging error message or removing the battery before the charging process has finished

To reduce the risk of damage to the plug and the power cord, disconnect the charger from the power supply by pulling on the plug and not on the cord.

The manufacturer recommends an annual service of the product (*applies to a device for which the expected lifetime has not yet been exceeded*).

3. CLEANING, DISINFECTION, STERILIZATION

3.1. GENERAL INFORMATION⁽²⁾

Power devices and attachments are often subject to high mechanical stresses and shocks during use and may be damaged. However, proper handling helps to prolong the life of surgical products.

ChM has validated the following instructions for reprocessing, fulfilling the requirements of EN ISO 17664. Observing the instructions provided below is important for the entire system to be clean. Nevertheless, other alternative methods of reprocessing can be equally effective.

Frequent reprocessing does not have a significant impact on the handpiece and attachments life. Gentle care can substantially increase the reliability and life of the system components.

The manufacturer shall not be liable for any damage resulting from unauthorized repairs of the product.

Cleaning and sterilization of products (*batteries and charger are not subject to sterilization*) **must take place before the first and every subsequent use.**

Reprocessing of the products must take place immediately after each use.

Cannulas, sleeves, collars, triggers and narrow places require special attention during cleaning.

Cleaning agents of pH 7 - 9.5 (neutral) which are commonly used in medicine for cleaning of medical devices (e.g. MediClean) are recommended. The use of cleaning agents with higher pH may, depending on the cleaning agent used, dissolve aluminum and its alloys, plastics or composite materials. PH values higher than 11 may also have a negative effect on stainless steel surfaces. Observe the instructions for use as well as information on the: concentration of diluted fluid, temperature, duration of exposure and water quality of enzymatic cleaning agents supplied by the manufacturers of these agents. If the temperature and time are not defined, follow the ChM recommendations. Products should be cleaned in freshly prepared solution.

Cleaning agents used will be in contact with the following materials: stainless steel, aluminum, plastics.

Personnel that re-processes the products should always be provided with personal protective garment in accordance with the information contained in the manual instructions and the safety data sheet of the product, supplied with the cleaning agent.

ChM recommends to use new sterile cutting elements every time.

⁽²⁾ Should the above-mentioned information be insufficient, please, contact the manufacturer.

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

For more information, please refer to the laws and guidelines for your country. Additional compliance with in-hospital policy and procedures, manufacturers recommendations of detergents, disinfectants and machines for clinical reprocessing are required.

The person that reprocesses the devices is responsible for ensuring that reprocessing performed has achieved the desired effect due to properly and correctly assembled, maintained, validated and approved equipment, as well as the use of appropriate materials, engagement of trained and authorized personnel, performed in the designated for reprocessing room. Any changes, deviations from the instructions introduced by the person responsible for the reprocessing process should be duly estimated in terms of their effectiveness and potential negative consequences.

3.2. HANDPIECE AND ATTACHMENTS

Prior to cleaning and disinfection, remove any attached accessories from the handpiece, such as attachments, cutting elements, battery, etc. All disposable instruments must not be reprocessed.

3.2.1. MANUAL CLEANING

Immediately after the surgery, the handpiece should be switched into a safety mode and the debris should be initially removed using disposable cleaning cloths. Do not use saline to moisten or soak the product before shipping to the decontamination room. Should the products be transported to the decontamination room with delay, they should be covered with a damp cloth. The wet cloth will prevent drying of the debris and will facilitate subsequent decontamination processes.

1. Rinse the product under running water for at least two minutes (*until visible dirt is removed*). Use a sponge, soft, lint-free cloth or soft bristled brush to assist in removing gross soil. Do not use for cleaning brushes with metallic bristles or pointed objects, as this may damage the outer coating of the product and consequently result in corrosion. Particular attention should be paid to cannulas, fixings, triggers and other gaps and attachments (*Fig.1, Fig.2, Fig.3*). Move all movable parts, such as triggers, locking sleeve for attachments to remove dirt.
2. Spray the device using a neutral pH enzymatic solution, leave for at least 2 minutes and wipe dry, thoroughly. The correct temperature, water quality and dilution/ concentration should be used in accordance with the enzymatic detergent manufacturer's directions.
3. Rinse the product under running water for at least 2 minutes, again. To clean cannulas and other hard-to-reach places. use a syringe, pipette, or water jet.
4. Again, the product should be cleaned under running water with an enzymatic neutral pH detergent (*e.g. MediClean*) for at least 5 min. For cleaning, use a sponge, a soft, lint-free cloth, or soft-bristle brushes. Move all moving parts, such as triggers, locking sleeve for attachments to remove dirt. Follow the manufacturer's instructions regarding the temperature, water quality and concentration/ dilution of the detergent. The correct temperature, water quality and dilution/ concentration should be used in accordance with the enzymatic detergent manufacturer's directions.
5. Rinse the product in lukewarm water for at least 2 minutes. To clean cannulas and other hard-to-reach places. use a syringe, pipette, or water jet. Move all moving parts, such as triggers, locking sleeve for attachments for thorough rinsing.
6. Visually inspect the device. Repeat steps 1–6 in freshly prepared enzymatic cleaning solution (*e.g. MediClean*)

until no visible dirt remains.

7. Finally, rinse the product in deionised water for at least 2 minutes.
8. The last step of manual cleaning is to dry the product with a soft, lint-free cloth or pure compressed medical air.
9. After the manual cleaning process, the functionality of the product should be tested.

Fig.1. Locking sleeve for attachments

Fig.2. Keyless attachment holder with max. open jaws which should be thoroughly cleaned

Fig.3. The battery cover should be removed during cleaning.

Do not immerse the handpiece in any liquid for longer than necessary. Do not soak the product in any liquid.

Do not wash the handpiece in the ultrasound cleaner - (does not apply to attachments).

Manual cleaning must always precede automatic cleaning/ disinfection.

3.2.2. WASHING - DISINFECTION IN A WASHER-DISINFECTOR

Having performed the manual cleaning (*initial*), the handpiece, attachments and cover must undergo the process of automatic washing - disinfection in a washer-disinfector. The washer-disinfector should meet the requirements of EN ISO 15883-2.

Other cleaning/disinfection procedures than the procedure described below (*including manual pre-cleaning*) **have not been validated by ChM.**

Place all devices in a specially designed tray for machine washing so that the elements do not come into contact with each other and do not move. Ensure that all cannulations (*handpiece and attachments*) are positioned vertically, i.e. in an upright position so that water outflow is facilitated.

Parameters of the automatic cleaning /disinfection cycle		
Phase	Duration (<i>minimum</i>)	Recommended
Initial rinsing	2 minutes	Cold water
Main washing	10 minutes	Warm water ($\geq 55^{\circ}\text{C}$) + MediClean (<i>pH</i> 7- 9.5)
Rinsing, Thermal disinfection	10 minutes	Deionized or demineralised hot water (90°C)
Drying	40 minutes	90°C

After the cleaning-disinfection is completed, the devices should be removed from the washer-disinfector and checked for dirt residue, especially in hard-to-reach areas such as cannulation, locking sleeve for attachment, etc. In case of dirt presence, manual cleaning /automatic cleaning-disinfection processes shall be repeated.

Prior to sterilization, devices must undergo visual inspection by trained personnel. Generally, unmagnified visual inspection under good light conditions is sufficient. All components must be inspected for dirt, damage and /or corrosion. Then the functional control of the device should be carried out, e.g. by operating the system in the air and checking all moving parts (*e.g. triggers, sleeves*) for smooth operation.

3.2.3. STERILIZATION

Packaging

The packaging process, prior to sterilization, should be validated by the end user (*health care facility*) according to requirements of EN ISO 11607-2 and should be carried out under controlled purity conditions. The washed, disinfected and dried product should be packaged in a steam sterilized packaging (*double foil-and-paper sleeve*) that complies with the requirements of EN ISO 11607-1 and is CE marked. The product must be packed in such a way that no contamination occurs at the time of unpacking the device for use. The sterilization packaging is intended to maintain the sterility of medical devices after the sterilization process during storage prior to use. Particular attention should be paid to the sharp and pointed ends of the instruments so as not to damage other adjacent components and not to break the sterile barrier of the packaging.

Sterilization

- For handpiece, attachments and cover sterilization, use the method provided below. Using another method may damage the components or result in ineffective sterilization of the products.

- The handpiece must be sterilized when open, i.e. the battery cover must be removed.
- Observe the recommended drying time to prevent accumulation of water inside the product. Moisture can cause corrosion and may interfere with proper sterilization.
- Always make sure that the products are completely dry before the sterilization process.
- **ChM** recommends sterilization using overheated steam in an autoclave. Nevertheless, the design and properties of a given autoclave can affect the end result of the sterilization process. Therefore, the healthcare facility should verify its current sterilization process using the current equipment and people who carry out the sterilization process, and validate it in accordance with requirements of EN ISO 17665-1.
- Responsibility for verification of the sterilization process is borne by the healthcare facility. To ensure the effectiveness of in-hospital sterilization of products, all cycles and methods should be verified using different sterilization chambers, products packaging and products distribution systems in the sterilization basket.
- The minimum drying time presented in the table below has been determined and validated for a single handpiece in the package. In case of a different sterilization configuration, e.g. several handpieces are packed or sterilization requires the use of a sterilization cassette or several palettes for instruments, drying time should increase, but the effectiveness of such sterilization should be verified.
- Drying time depends on factors such as: product distribution, tray material, packaging method. For the handpiece and attachments, the minimum drying time is 20 minutes for each component individually.
- Once the sterilization process has been completed, the products should "cool down".

Validated parameters of overheated steam sterilization cycles:

- temperature: 134°C
- minimum exposure time: 5 min
- minimum drying time: 20 min

**The following maximum values may not be exceeded:
140°C over a maximum of 22 minutes.**

Higher values can damage the sterilized products.

- Hot air sterilization, ethylene oxide and formaldehyde sterilization are not recommended.
- Handpiece reliability is guaranteed for 100 full re-sterilization cycles (*including cleaning, disinfection and sterilization*).

3.3. STORAGE AND HANDLING BEFORE USE

- Always transport the packaged articles with great care to avoid damaging the packaging sterile barrier.
- Packaged, sterilized products should be stored in a controlled environment where they will not be exposed to high /low temperature and humidity.
- Excessive displacement of sterilized, packaged items can breach their sterile barrier and lead to contamination.
- **CharDrive** system should be stored in appropriate protective packaging, in a clean, dry place, at room temperature and under conditions that provide protection from direct sunlight.

The device must not be stored or operated in explosive atmospheres.

3.4. BATTERIES AND BATTERY CHARGER

The battery immediately after the surgery should be removed from the handpiece. Verify it for cracks or damage. Do not use the battery if damaged.

The battery and charger must not be sterilized, washed, rinsed, soaked or dropped onto the ground. These processes can destroy the above items as well as cause secondary damage.

The battery can only be cleaned with a clean, dry cloth, without any detergent.

After each use, the battery should be placed in the charger for recharging.

The charger should be disconnected from the power supply, and then its surface may be wiped with a soft cloth moistened but not dripping with mild disinfectant and dried immediately.

Do not allow water to come into contact with the charger's charging ports as this could lead to product damage.

4. REPAIRS

The device should be sent to **ChM** for repair if malfunctions or recognized as damaged, to the following address:

ChM sp. z o.o.
Lewickie 3b
16-061 Juchnowiec Kościelny, POLAND
tel: +48 85 86 86 100
fax: +48 85 86 86 101
e-mail: chm@chm.eu

- If the product has been dropped on a hard surface, it should be sent to the manufacturer.
- Damaged products must not be used. If repair is not possible, the product should be disposed of in accordance with the DISPOSAL section.

The manufacturer assumes no liability for any damage resulting from unauthorized repairs of the product.

5. LIFETIME

The lifetime of the **CharDrive** Handpiece is assumed for 10 years. After this period, the manufacturer (**ChM** sp. z o.o.) may refuse to perform service repairs or technical inspections. The battery is a consumable component. Its life is determined by the way and frequency of use. The manufacturer recommends replacing batteries with new ones at least once a year. Handpiece reliability is guaranteed by 100 full re-sterilization cycles (*including cleaning, disinfection and sterilization*).

6. DISPOSAL

In most cases, faulty devices can be repaired (*refer to chapter REPAIRS*).

Rechargeable lithium-ion battery.

Under Directive 2006/66/EC, batteries must be collected and brought into use in accordance with the rules in force in the EU Member States.

Risk of fire, explosion and burns. DO NOT disassemble, crush, heat above 100 °C, or incinerate the battery cells.

Contaminated products must be full reprocessed so that the risk of infection during product disposal will be excluded.

The product may not be disposed of with household waste.

7. TROUBLESHOOTING ⁽³⁾

7.1. CharDrive HANDPIECE

Problem	Possible causes	Solution
Handpiece does not start.	No battery in handpiece	Insert charged battery.
	The battery is flat	Charge the battery
	Safety mode is activated	Set the selected mode of operation
	Thermal overload protection system is activated	Wait until cooled down
	Handpiece is faulty	Send to ChM
Handpiece does not have enough power	The battery is flat	Charge the battery
	Handpiece is faulty	Send to ChM
Handpiece stops suddenly	The battery is flat	Charge the battery
	Thermal overload protection system is activated	Wait until cooled down
	Handpiece is faulty	Send to ChM
Handpiece continues to operate after releasing the trigger	Trigger is blocked by surgery debris	Press trigger several times; clean as recommended
	Handpiece is faulty	Send to ChM
Handpiece heats up.	Overloaded device	Wait until cooled down
	Non-observance of duty cycles of the device	Observe the duty cycles of the device
	Handpiece is faulty	Send to ChM
Attachments cannot be attached/ detached to/ from the handpiece.	Plugged locking sleeve for attachments	Clean as recommended
	Handpiece is faulty	Send to ChM
	Attachment is faulty	Send to ChM
	Plugged attachment	Clean as recommended
The triggers are difficult to move	Plugged triggers	Press trigger several times; clean as recommended
		Turn off all electrical equipment not in use
Handpiece self-activates	Operating the handpiece near a device that emits a strong electromagnetic field	Change the location of electrical equipment
		Plug in electrical equipment to other available outlets
Inadequate direction of rotation of the handpiece	Device control program failure	Send to ChM
Handpiece overheated	Negligence of Instructions for Use supplied with the product	Consult Instructions for Use of the product
	Non-observance of duty cycles of the device	Observe the duty cycles of the device
Handpiece run-out, in extreme cases, handpiece being torn from the operator's hand	Improper installation of the cutting element, e.g. a drill in the attachment	Re-install the cutting element in the attachment
	Improper installation of the attachment to the handpiece	Re-install the attachment in the handpiece
	Handpiece is faulty	Send to ChM

(3) If any other problems than the ones listed below have occurred or if the suggested actions have not solved the problem, please contact the manufacturer. DO NOT perform repairs on your own unless otherwise stated in this manual.

Problem	Possible causes	Solution
Battery cover cannot be fitted onto the handpiece.	Insufficient / incorrect placement of the battery in the handpiece	Press tight/ place the battery properly in the handpiece
	Use of non-original battery	Use only original batteries
No functionality of DRILL and / or REAM modes	Inadequate re-processing of the handpiece which can result in dirt drying off	Clean as recommended
	Handpiece is faulty	Send to ChM
No functionality of the LED system informing about the battery charge level in the handpiece	Incorrect re-processing of the handpiece	Contact the manufacturer
	No battery in the handpiece	Place the charged battery in the handpiece
	Flat battery in the handpiece	Place the charged battery in the handpiece
	Faulty LED system	Contact the manufacturer

7.2. CharDrive BATTERY

Problem	Possible causes	Solution
Battery cover cannot be fitted onto the handpiece.	Insufficient / incorrect placement of the battery in the handpiece	Press tight/ place the battery properly in the handpiece
	Use of non-original battery	Use only original batteries
The battery does not fit into the charger port	Use of non-original battery	Use only original batteries
Excessive heating up of the battery during the charging process	<i>Battery - charger incompatibility</i>	Use only original batteries/ charger
	<i>Damaged battery and/ or charger (e.g. by washing, sterilization)</i>	Send to ChM
Battery got discharged while performing the surgery	<i>Excessive use of the battery</i>	Recharge the battery
	Not fully charged battery	Recharge the battery
	Battery is faulty	Send to ChM
	Battery has been re-processed	Send to ChM
Fully charged battery does not work	Battery has been re-processed	Send to ChM
	Battery has been dropped onto a hard surface	Send to ChM
Accelerated battery degradation	Working in conditions harmful to the battery	Send to ChM
	Deep discharge of the battery	Send to ChM
	Inadequate storage and in-hospital transportation	Send to ChM
The battery and / or charger has been washed/ sterilized	Staff negligence	Send to ChM
	Negligence of Instructions for Use supplied with the product	Send to ChM

7.3. ATTACHMENTS AND CUTTING ELEMENTS

Problem	Possible causes	Solution
Attachments cannot be attached/ detached to/ from the handpiece.	Plugged locking sleeve for attachments	Clean as recommended
	Handpiece is faulty	Send to ChM
	Attachment is faulty	Send to ChM
	Plugged attachment	Clean as recommended

Problem	Possible causes	Solution
Handpiece run-out, in extreme cases, handpiece being torn from the operator's hand	Improper installation of the cutting element, e.g. a drill in the attachment	Re-install the cutting element in the attachment
	Improper installation of the attachment to the handpiece	Re-install the attachment in the handpiece
	Handpiece is faulty	Send to ChM
The handpiece attachment/ cutting element does not rotate despite the activated motor	Handpiece is faulty	Send to ChM
	Improper installation of the attachment to the handpiece	Re-install the attachment in the handpiece
	Applying excessive pressure	Apply less pressure
Excessive increase in tissue and/ or cutting element temperature	Blunt cutting element	Replace the instrument for new one
	Non-observance of duty cycles of the device	Observe the duty cycles of the device
Cutting element cannot be installed in the attachment	Worn, damaged attachment	Send to ChM
	Worn cutting element	Replace the instrument for new one
	Incompatible cutting element tip and attachment	Choose a cutting element that is compatible with the attachment or another way round

7.4. CharDrive BATTERY CHARGER

Problem	possible causes	Solution
Electrical interference	Occurrence of electrical noise	Turn off all electrical equipment not in use
		Change the location of electrical equipment
		Plug in electrical equipment to other available outlets
The battery does not fit into the charger port	Use of non-original battery	Use only original batteries
Excessive heating up of the battery during the charging process	<i>Battery - charger incompatibility</i>	Use only original batteries/ charger
	<i>Damaged battery and/ or charger (e.g. by washing, sterilization)</i>	Send to ChM
The battery and / or charger has been washed/ sterilized	Staff negligence	Send to ChM
	Negligence of Instructions for Use supplied with the product	Send to ChM
No messages displayed by the charger	The charger is not connected to the power source	Connect the charger to a power source
	Damage to the charger (e.g. as a result of dropping, reprocessing)	Send to ChM
A warning appears when charging the battery	Damaged battery /charger	Send to ChM
	Battery has been removed from the charger before its 100% charged	Charge the battery fully
Loosened the charger port used to charge the battery	Loosened the screw securing the battery charging socket	Tighten the screw
Possible electrical shock	Damage to the power cord of the charger	Send to ChM

8. TECHNICAL DATA ⁽⁴⁾

8.1. DUTY CYCLE

Intermittent operation type S9, according to EN 60034-1 

X_{on}	Y_{off}	Cycles
30 sec.	60 sec.	5

Electrical systems may heat up if in constant use. Therefore, the handpiece should be allowed to cool for at least 60 seconds following 30 seconds of constant use. After 5 such cycles, the handpiece should be allowed to cool for 30 minutes. If the above-mentioned is observed, the system will be prevented from overheating and the patient or user from possible injury. The user is responsible for the use of the system as indicated. If longer periods of constant use are required, an additional handpiece should be available.

- **Observe the prescribed duty cycles, carefully.**
- **Always use new cutting instruments to prevent heating up of the system due to reduced cutting characteristics.**
- **Cool manually instruments with irrigation liquid to prevent heat necrosis of tissues.**
- **The battery charger is designed for continuous operation.**

⁽⁴⁾ Technical data is subject to tolerances.

8.2. CharDrive HANDPIECE

Dimensions of handpiece without attachment	223x41x220mm
Weight of handpiece	1330g
Gradually adjustable speed	mode DRILL: 0–1000 rpm mode REAM: 0–250 rpm
Handpiece cannulation	3.4mm
Protection class	BF, EN 60601-1
Degree of protection	IP46, EN 60529
Directive 93/42/EEC classification	Ila
Power supply	Internally powered with lithium-ion batteries
Number of triggers	2
EMC	EN 60601-1-2
Applicable standards	EN 60601-1

8.3. CharDrive BATTERY

Type	Li-Ion
Dimensions	142x42x26mm
Weight	205g
<i>Operating voltage (normal)</i>	11.1 VDC
Capacity	2.2 Ah
Typical charging time	≤ 90 min
Directive 93/42/EEC classification	I
EMC	EN 60601-1-2
Applicable standards	EN 60601-1

8.4. CharDrive BATTERY CHARGER

Maximum input current	0.65A
Power cord	CEE 7/7 (E/F) angle plug - IEC C13
<i>Dimensions (LxWxH)</i>	136x244x127mm
Weight	2130g
EMC	EN 60601-1-2
Applicable standards	EN 60601-1
Degree of protection	IP30, EN 60529
Directive 93/42/EEC classification	I
Protection class	I class

9. ELECTROMAGNETIC COMPATIBILITY IN ACCORDANCE WITH EN 60601-1-2

The **CharDrive** Medical Drive System (*handpiece and battery charger*) is intended for use in the electromagnetic environment specified below. The customer or the user of the System should assure that it is used in such an environment.

9.1. CharDrive HANDPIECE

9.1.1. ELECTROMAGNETIC EMISSIONS

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions, CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emission is very low and it is not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions, IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions, IEC 61000-3-3	Not applicable	

9.1.2. ELECTROMAGNETIC IMMUNITY

Immunity test standard	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±8 kV contact ±15 kV air	Floors should be made from wood, concrete or ceramic tile. If floors are covered with a synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for in/ out lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% UT (0.5 cycle) 40% UT (5 cycles) 70% UT (25 cycles) <5% UT for 5 s	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires its continuous use even during power outages, it is recommended that the device be connected to an emergency power supply.
Power frequency (50/60 Hz) magnetic field (50/60Hz) IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Immunity test standard	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 800 MHz	Nie dotyczy	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 800 MHz	E1 = 10 V/m	Recommended separation distance $d = 1.17 \sqrt{P}$ $d = 0.35 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 0.70 \sqrt{P} \text{ 800 MHz to 2.5 GHz}$ where P is the maximum output power rating of the transmitter in watts (<i>W</i>) according to the transmitter manufacturer and d is the recommended separation distance in metres (<i>m</i>). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
	3 V/m 80 MHz to 800 MHz	E2 = 10V/m	

- a) Field strengths from fixed transmitters, such as base stations for radio telephones (*cellular/cordless*) and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

9.1.3. SEPARATION DISTANCES

Rated maximum output power of transmitter (<i>W</i>)	Separation distance according to frequency of transmitter (<i>m</i>)		
	150 kHz to 80 MHz $d = 1.17 \sqrt{P[m]}$	80 MHz to 800 MHz $d = 1.17 \sqrt{P[m]}$	800 MHz to 2.5 GHz $d = 2.33 \sqrt{P[m]}$
0.01	0.12	0.04	0.07
0.1	0.37	0.11	0.22
1	1.2	0.35	0.7
10	3.7	1.1	2.2
100	12	3.5	7

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in metres (*m*) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (*W*) according to the transmitter manufacturer.

9.2. CharDrive BATTERY CHARGER

9.2.1. ELECTROMAGNETIC EMISSIONS

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions, CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emission is very low and it is not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions, IEC 61000-3-2	Class A	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations / flicker emissions, IEC 61000-3-3	Complies	

9.2.2. ELECTROMAGNETIC IMMUNITY

Immunity test standard	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be made from wood, concrete or ceramic tile. If floors are covered with a synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for in/ out lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	±1 kV line to line ±2 kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% UT (0.5 cycle) 40% UT (5 cycles) 70% UT (25 cycles) <5% UT for 5 s	<5% UT (0.5 cycle) 40% UT (5 cycles) 70% UT (25 cycles) <5% UT for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires its continuous use even during power outages, it is recommended that the device be connected to an emergency power supply.
Power frequency (50/60 Hz) magnetic field (50/60Hz) IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Immunity test standard	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 800 MHz	V1 = 10Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 800 MHz	E1 = 10 V/m	Recommended separation distance $d = 1.17 \sqrt{P}$ $d = 0.35 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 0.70 \sqrt{P} \text{ 800 MHz to 2.5 GHz}$ where P is the maximum output power rating of the transmitter in watts (<i>W</i>) according to the transmitter manufacturer and d is the recommended separation distance in metres (<i>m</i>). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
	3 V/m 80 MHz to 800 MHz	E2 = 7V/m	

- a) Field strengths from fixed transmitters, such as base stations for radio telephones (*cellular/cordless*) and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

9.2.3. SEPARATION DISTANCES

Rated maximum output power of transmitter (<i>W</i>)	Separation distance according to frequency of transmitter (<i>m</i>)		
	150 kHz to 80 MHz $d = 1.17 \sqrt{P[m]}$	80 MHz to 800 MHz $d = 1.17 \sqrt{P[m]}$	800 MHz to 2.5 GHz $d = 2.33 \sqrt{P[m]}$
0.01	0.12	0.04	0.1
0.1	0.37	0.11	0.32
1	1.2	0.35	1.0
10	3.7	1.1	3.2
100	12	3.5	10

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in metres (*m*) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (*W*) according to the transmitter manufacturer.

10. EXPLANATION OF SYMBOLS USED

Symbol	Meaning	Description
	Do not discard the product with household waste	Worn devices and other electrical products must be collected separately and disposed of in accordance with the applicable electrical and electronic equipment law.
	Keep this product dry!	-
	Classification Type BF	Degree of protection against electric shock.
	European conformity mark with the notified body number	Conform to the European Medical Device Directive 93/42/CEE for class IIa medical devices (<i>notified body number</i>)
	European conformity mark	Conform to the European Medical Device Directive 93/42/CEE for class I medical devices
	Instructions for Use	Prior to first use, consult Instructions for Use
IPxx	Protection degree	Ingress protection degree
	Caution	Read the provided Instructions for Use before operating the device.
	Manufacturer/ Date of manufacture	The address of the device manufacturer and the date of manufacture of the device.

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