


ChM[®]

7.0 ChM Locked Plating
ChLP system


































7.0ChLP femoral periprosthetic plate

3.7220; 3.7221;
3.7222; 3.7223;
3.7224; 3.7225;
3.7276; 3.7277

- *IMPLANTS*
- *INSTRUMENT SET*
- *SURGICAL TECHNIQUE*



SYMBOLS DESCRIPTION

	Titanium or titanium alloy		H length [mm]
	Cobalt		Angle
	Left		available lengths
	Right		Available number of holes
	Available versions: left/right		Thickness [mm]
	Length		Scale 1:1
	Torx drive		Number of threaded holes in the shaft part of the plate
	Torx drive cannulated		Number of locking holes in the plate
	Hexagonal drive		Variable angle
	Hexagonal drive cannulated		Cortical
	Cannulated		Cancellous
	Locking		Available in sterile/ non- sterile condition
	Diameter [mm]		Refer to surgical technique
	Caution - pay attention to a special procedure.		
	Perform the activity under X-Ray control.		
	Information about the next stages of a procedure.		
	Proceed to the next stage.		
	Return to the specified stage and repeat the activity.		
	Before using the product, carefully read the Instructions for Use. It contains, among others, indications, contraindications, side effects, recommendations and warnings related to the use of the product.		
	The above description is not a detailed instruction of conduct. The surgeon decides about choosing the operating procedure.		

www.chm.eu

Document No ST/80-703
 Date of issue 12.06.2019
 Review date P-004-30.08.2023

The manufacturer reserves the right to introduce design changes.

Updated INSTRUCTIONS FOR USE are available at the following website: ifu.chm.eu

1. INTRODUCTION	5
2. IMPLANT DESCRIPTION	6
3. SURGICAL TECHNIQUE - 7.0ChLP FEMORAL PERIPROSTHETIC PLATE	9
3.1. PATIENT'S POSITIONING	9
3.2. SURGICAL APPROACH	9
3.3. FRACTURE REDUCTION	9
3.4. Implant selection	9
3.5. PLATE INSERTION	9
3.6. TEMPORARY PLATE STABILIZATION	10
3.7. SCREWS INSERTION	10
3.8. WOUND CLOSURE	10
4. SURGICAL PROCEDURES	11
4a. PROCEDURE OF TEMPORARY IMPLANT STABILIZATION	11
4b. PROCEDURE OF CORTICAL SELF-TAPPING SCREW 4.5 [3.1471] INSERTION	12
4c. PROCEDURE OF 7.0ChLP SELF-TAPPING SCREW 5.0 [3.5210] INSERTION	14
4d. PROCEDURE OF 7.0ChLP PERIPROSTHETIC SCREW 5.0 [3.5247] INSERTION	16
4e. PROCEDURE OF 7.0ChLP SCREW VA 4.0 [4.5246] INSERTION	17
4f. PROCEDURE OF 7.0ChLP SCREW VA 5.0 [4.5244] INSERTION	19
5. POSTOPERATIVE PROCEDURE	21
6. IMPLANT REMOVAL	21
7. CATALOGUE PAGES	22
7a. INSTRUMENT SET	22
7b. Implants	24
7c. Screws	29

1. introduction

This surgical technique applies to 7.0ChLP locked plating system used for stabilization of proximal femur fractures. The plates are a part of the ChLP locked plating system developed by **ChM**. The presented range of implants is made of materials in accordance with ISO 5832 standards.

The system includes:

- implants (*plates and screws*),
- instrument set used in the surgery,
- surgical technique.

Indications

- trochanteric osteotomies,
- trochanteric fractures,
- periprosthetic femur fractures.

Plate selection and shaping

The plates are available in various lengths and for left and right limb separately. This allows for optimal selection of the implant to the fracture type. Shaping of the plates in their epiphyseal part is not allowed.




Before using the product, carefully read the Instructions for Use. It contains, among others, indications, contraindications, side effects, recommendations and warnings related to the use of the product.




The above description is not a detailed instruction of conduct. The surgeon decides about choosing the operating procedure.


2. IMPLANT DESCRIPTION

Femoral periprosthetic plates are a part of 7.0ChLP system. This system includes also compatible locking screws. To facilitate their identification, both titanium plate and screws are blue anodized.



7.0ChLP femoral periprosthetic plate






3. SURGICAL TECHNIQUE - 7.0ChLP femoral periprosthetic plate


3.1. PATIENT'S POSITIONING

Place a patient supine. Ensure lateral and AP visualization.



3.2. SURGICAL APPROACH

Lateral access. Perform a more or less extensive incision of the skin (*depending on the implant used*). The incision shall start from the top of the greater trochanter to the lateral condyle of the femur.



3.3. FRACTURE REDUCTION


Perform fracture reduction. If need be, temporarily stabilize the bone fragments with Kirschner wires and/or reduction pliers.

3.4. Implant selection

Select the right size of an implant to the type of fracture, bone size and structure.

3.5. PLATE INSERTION

Position the implant correctly on the bone.





3.6. TEMPORARY PLATE STABILIZATION

Stabilize the position of the implant inserting Kirschner wires into appropriate holes (acc. to procedure 4a).




The Kirschner wires in the epiphyseal part indicate (in the transverse view) the plane of the screws to be inserted. They facilitate positioning of the plate over the prosthesis.



3.7. SCREWS INSERTION

Insert screws of a suitable length into the holes of the plate.

- Cortical self-tapping screw 4.5 [3.1471] (acc. to procedure 4b).
- 7.0ChLP self-tapping screw 5.0 [3.5210] (acc. to procedure 4c).
- 7.0ChLP periprosthetic screw 5,0 [3.5247] (acc. to procedure 4d).
- 7.0ChLP screw VA 4.0 [4.5246] (acc. to procedure 4e).
- 7.0ChLP screw VA 5.0 [4.5244] (acc. to procedure 4f).



Insert the cortical screws 4.5 into a bone fragment before inserting the locking screws.



The doctor decides about the order and number of locking and cortical screws to be inserted.

3.8. WOUND CLOSURE


Before closing the wound, take an X-Ray image in at least two projections to confirm implant position and fracture reduction. Make sure all the screws are properly tightened and do not penetrate the joint surface. Use appropriate surgical technique to close the wound.

4. SURGICAL PROCEDURES

4a. PROCEDURE OF TEMPORARY IMPLANT STABILIZATION



Stabilization using Kirschner wires

- Stabilize temporary the implant inserting Kirschner wires 2.0/210 [40.4815.210] into dedicated holes in the plate.





Stabilization in locking holes using Kirschner wires

- Insert guide sleeve 7.0/4.0 [40.5705.740] into the locking hole of the plate.
- Insert Kirschner wire [40.4815.210] through the guide sleeve 7.0/4.0 [40.5705.740].



Stabilization using setting-compressing screw


- Insert guide sleeve 7.0/4.0 [40.5705.740] into the locking hole of the plate.
- Insert setting-compressing screw 4.0/180 [40.5706.740] through the guide sleeve 7.0/4.0 [40.5705.740].
- Tighten the nut of the setting-compressing screw 4.0/180 [40.5706.740] and push the plate to the bone.



4b. PROCEDURE OF CORTICAL SELF-TAPPING SCREW 4.5 [3.1471] INSERTION


Compression guide positioning

Position the compression guide VA 4.0 [40.8207.040] in a desired position:




Hole drilling

Perform a hole through both cortices for a cortical screw 4,5 insertion. For drilling, use drill with scale 3,2/210 [40.5650.212] and compression guide in a desired position.




Measurement of hole depth

OPTION 1: Determine the length of the screw using the scale on the drill with scale 3,2/210 [40.5650.212].





OPTION II: Insert depth measure [40.4639.550] into drilled hole until the hook of the measure rests against the outer surface of the second cortex.



Screw insertion


Insert cortical screw using torque limiting ratchet T handle 4Nm [40.6660.000] and screwdriver tip T25 [40.5684.200].



4c. PROCEDURE OF 7.0ChLP SELF-TAPPING SCREW 5.0 [3.5210] INSERTION

Guide sleeve insertion

- Insert guide sleeve 7.0/4.0 [40.5705.740] into the locking hole of the plate.




Hole drilling

Drill using drill with a scale 4.0/210 [40.5651.212] until desired depth is reached.





Measurement of hole depth



OPTION I: Determine the length of the screw using the scale on the drill with scale 4.0/210 [40.5651.212].



OPTION II: or using measure on VA 4.0 guide.





OPTION III: Having removed the guide sleeve 7.0/4.0 [40.5705.740], use depth measure [40.4639.550] to determine the length of the screw.



Screw insertion


Remove guide sleeve 7.0/4.0 [40.5705.740]. Insert locking screw using torque limiting ratchet T handle 4Nm [40.6660.000] and screwdriver tip T25 [40.5684.000].



4d. PROCEDURE OF 7.0ChLP PERIPROSTHETIC SCREW 5.0 [3.5247] INSERTION

Guide sleeve insertion

- Insert guide sleeve 7.0/4.0 [40.5705.740] into the locking hole of the plate.



Hole drilling


Drill using drill with a scale 4.0/210 [40.5651.212] until desired depth is reached.




Drill under X-Ray control to avoid a collision of the drill with already implanted screws.

Measurement of hole depth

OPTION I: Determine the length of the screw using the scale on the drill with scale 4.0/210 [40.5651.212].





OPTION II: Having removed the guide sleeve 7.0/4.0 [40.5705.740], use depth measure [40.4639.550] to determine the length of the screw.



Screw insertion

Remove guide sleeve 7.0/4.0 [40.5705.740]. Insert locking screw using torque limiting ratchet T handle 4Nm [40.6660.000] and screwdriver tip T25 [40.5684.000].



4e. PROCEDURE OF 7.0ChLP SCREW VA 4.0 [4.5246] INSERTION



When using variable angle (VA) screws, there is a risk of collision of screws or a drill with already implanted screws. Well-thought-out trajectory of inserted screws and intraoperative X-Ray control of drilling reduces the risk of the collision.

Guide VA positioning

- Insert the guide VA 4,0 [40.8207.040] na pełną głębokość w osi otworu blokowanego. into the locking hole co-axially.
- Set the desired inclination of the guide in relation to the locking hole axis. The guide enables the inclination of 15° in each direction with respect to the axis of the locking hole.



40.8207.040



Exceeding the inclination angle of more than 15° may prevent proper locking of the VA screw in the plate hole.

Hole drilling

- Drill using drill with scale 3,2/210 [40.5650.212] until desired depth is reached.



40.5650.212




Drill under X-Ray control to avoid a collision of the drill with already implanted screws.

Measurement of hole depth


OPTION I: Read the length of the screw from the drill measure [40.5650.212]



40.5650.212





OPTION II: Having removed the guide VA, use depth measure **[40.4639.550]** to determine the length of the screw.



Screw insertion

Use torque limiting ratchet handle 4Nm **[40.6660.000]** and screwdriver tip T25 **[40.5684.200]** to insert the VA screw.



Change of the VA screw positioning

It is possible to lock the VA screw three times in the threaded hole of the plate.
The hole in the plate in which the VA screw was locked cannot be used to insert a standard locking screw.

4f. PROCEDURE OF 7.0ChLP SCREW VA 5.0 [4.5244] INSERTION



When using variable angle (VA) screws, there is a risk of collision of screws or a drill with already implanted screws. Well-thought-out trajectory of inserted screws and intraoperative X-Ray control of drilling reduces the risk of the collision.

Guide VA positioning

- Insert the guide VA 4,0 [40.8207.040] into the locking hole co-axially.
- Set the desired inclination of the guide in relation to the locking hole axis. The guide enables the inclination of 15° in each direction with respect to the axis of the locking hole.



40.8207.040



Exceeding the inclination angle of more than 15° may prevent proper locking of the VA screw in the plate hole.

Hole drilling

- Drill using drill with scale 4,0/210 [40.5651.212] until desired depth is reached.



40.5651.212




Drill under X-Ray control to avoid a collision of the drill with already implanted screws.

Measurement of hole depth


WARIANT I: Read the length of the screw from the drill measure 4,0/210 [40.5651.212]



40.5651.212





WARIANT II: Having removed the guide VA, use depth measure **[40.4639.550]** to determine the length of the screw.



Screw insertion

Use torque limiting ratchet handle 4Nm **[40.6660.000]** and screwdriver tip T25 **[40.5684.200]** to insert the VA screw.



Change of the VA screw positioning

It is possible to lock the VA screw three times in the threaded hole of the plate.
The hole in the plate in which the VA screw was locked cannot be used to insert a standard locking screw.

5. POSTOPERATIVE PROCEDURE

Introduce appropriate postoperative treatment. The physician decides on the post-operative treatment and its conduct. In order to avoid patient's movement limitations, introduce exercises as soon after surgery as possible. However, make sure that the limb is not fully loaded before fragments osteosynthesis is complete.

6. IMPLANT REMOVAL


















The physician decides about implant removal. In order to remove the implants from the body, unlock all the locking screws first and then remove them from the bone. This will prevent any rotation of the plate when removing the last locking screw.

7. CATALOGUE PAGES

7a. INSTRUMENT SET

Instrument set for 7.0ChLP 4x4 1/2H

15.0207.208

	Name	Catalogue no.	Pcs.
	Tray for 7.0ChLP instrument set 4x4 1/2H	14.0207.208	1
	Kirschner wire 2.0/210	40.4815.210	4
	Drill with scale 3.2/210	40.5650.212	2
	Drill with scale 4.0/210	40.5651.212	2
	Cannulated drill with scale 5.0/2.2/210	40.5652.212	1
	Setting-compressing screw 4.0/180	40.5706.740	1
	Guide VA 4.0	40.8207.040	1
	Guide sleeve 7.0/4.0	40.5705.740	3
	Guide sleeve 7.0/3.2	40.5705.732	2
	Guide sleeve 9/5.0	40.5654.750	1
	Guide sleeve 5.0/2.0	40.5654.120	1
	Protective guide 9/7	40.5708.000	2
	Torque limiting ratchet handle T 4Nm	40.6660.000	1
	Screwdriver tip T25-1/4	40.5684.200	1
	Cannulated screwdriver tip T30-1/4	40.5685.200	1
	Depth measure	40.4639.550	1
Optional instrument			
	Torque connector 4Nm	40.5927.040	

OPTIONAL TOOLS


Optional instrument 7.0ChLP cerclage screws 3.1221.170



Tripod screwdriver tip 7.0ChLP


40.6271.700







7.0ChLP femoral periprosthetic plate



①		Ti		SHORT
Len	L	R		
6	222	3.7221.606	3.7220.606	
8	274	3.7221.608	3.7220.608	
10	326	3.7221.610	3.7220.610	
12	378	3.7221.612	3.7220.612	

②		LONG	
Len	L	R	
6	229	3.7223.606	3.7222.606
8	281	3.7223.608	3.7222.608
10	333	3.7223.610	3.7222.610
12	384	3.7223.612	3.7222.612



	TIA	Co	⊗	⊙	VA	↻	↗
 3.5210.xxx			✓	✓		✓	5.0
 3.5247.xxx			✓	✓		✓	5.0
 3.1471.xxx			✓			✓	4.5
 4.5244.xxx		✓	✓	✓	✓	✓	5.0
 4.5246.xxx		✓	✓	✓	✓	✓	4.0
 3.1221.170			✓				

7.0ChLP trochanteric periprosthetic plate



①		Ti	SHORT
Len	LR		
4	169	3.7276.604	
②			LONG
4	176	3.7277.604	

	TiA	Co	⊗	⊙	VA	C	↗
	3.5210.xxx		✓	✓		✓	5.0
	3.5247.xxx		✓	✓		✓	5.0
	3.1471.xxx		✓			✓	4.5
		4.5244.xxx	✓	✓	✓	✓	5.0
		4.5246.xxx	✓	✓	✓	✓	4.0
	3.1221.170			✓			



7.0ChLP trochanteric periprosthetic plate



①			Ti	SHORT
	Len		L R	
2	42		3.7224.600	

②				LONG
2	50		3.7225.600	



	TiA	Co			VA		
	3.5210.xxx		✓	✓		✓	5.0
	3.5247.xxx		✓	✓		✓	5.0
	3.1471.xxx		✓			✓	4.5
		4.5244.xxx	✓	✓	✓	✓	5.0
		4.5246.xxx	✓	✓	✓	✓	4.0
	3.1221.170			✓			

7.0ChLP femoral periprosthetic plate




	Len	L	Ti	R
6	196	3.7273.606		3.7274.606
8	248	3.7273.608		3.7274.608
10	300	3.7273.610		3.7274.610
12	352	3.7273.612		3.7274.612



	TiA	Co			VA		
	3.5210.xxx		✓	✓		✓	5.0
	3.5247.xxx		✓	✓		✓	5.0
	3.1471.xxx		✓			✓	4.5
		4.5244.xxx	✓	✓	✓	✓	5.0
		4.5246.xxx	✓	✓	✓	✓	4.0
	3.1221.170			✓			

7.0ChLP femoral periprosthetic plate



4	143	3.7275.604

	3.5210.xxx		✓	✓		✓	5.0
	3.5247.xxx		✓	✓		✓	5.0
	3.1471.xxx		✓			✓	4.5
		4.5244.xxx	✓	✓	✓	✓	5.0
		4.5246.xxx	✓	✓	✓	✓	4.0
	3.1221.170			✓			



7.0ChLP self-tapping screw 5.0



Len	TiA
16	3.5210.016
18	3.5210.018
20	3.5210.020
22	3.5210.022
24	3.5210.024
26	3.5210.026
28	3.5210.028
30	3.5210.030
32	3.5210.032
34	3.5210.034
36	3.5210.036
38	3.5210.038
40	3.5210.040
42	3.5210.042
44	3.5210.044
46	3.5210.046
48	3.5210.048
50	3.5210.050
52	3.5210.052
54	3.5210.054
56	3.5210.056
58	3.5210.058
60	3.5210.060
65	3.5210.065
70	3.5210.070
75	3.5210.075
80	3.5210.080
85	3.5210.085
90	3.5210.090
95	3.5210.095
100	3.5210.100
105	3.5210.105
110	3.5210.110

Cortical self-tapping screw 4.5



Len	TiA
16	3.1471.016
18	3.1471.018
20	3.1471.020
22	3.1471.022
24	3.1471.024
26	3.1471.026
28	3.1471.028
30	3.1471.030
32	3.1471.032
34	3.1471.034
36	3.1471.036
38	3.1471.038
40	3.1471.040
42	3.1471.042
44	3.1471.044
46	3.1471.046
48	3.1471.048
50	3.1471.050
52	3.1471.052
54	3.1471.054
56	3.1471.056
58	3.1471.058
60	3.1471.060
65	3.1471.065
70	3.1471.070
75	3.1471.075
80	3.1471.080
85	3.1471.085
90	3.1471.090
95	3.1471.095
100	3.1471.100
105	3.1471.105
110	3.1471.110

7.0ChLP periprosthetic screw 5,0



Len	TiA
10	3.5247.010
12	3.5247.012
14	3.5247.014
16	3.5247.016
18	3.5247.018
20	3.5247.020

7.0ChLP cerclage screws



TiA
3.1221.170



NOTE: Screw dedicated to periprosthetic plates.
Optional instruments for cerclage screws - Tripod screwdriver tip 7.0ChLP [40.6271.700]

7.0ChLP screw VA 4.0



Len	Co
16	4.5246.016
18	4.5246.018
20	4.5246.020
22	4.5246.022
24	4.5246.024
26	4.5246.026
28	4.5246.028
30	4.5246.030
32	4.5246.032
34	4.5246.034
36	4.5246.036
38	4.5246.038
40	4.5246.040
42	4.5246.042
44	4.5246.044
46	4.5246.046
48	4.5246.048
50	4.5246.050
52	4.5246.052
54	4.5246.054
56	4.5246.056
58	4.5246.058
60	4.5246.060
65	4.5246.065
70	4.5246.070
75	4.5246.075
80	4.5246.080
85	4.5246.085
90	4.5246.090
95	4.5246.095

7.0ChLP screw VA 5.0



Len	Co
16	4.5244.016
18	4.5244.018
20	4.5244.020
22	4.5244.022
24	4.5244.024
26	4.5244.026
28	4.5244.028
30	4.5244.030
32	4.5244.032
34	4.5244.034
36	4.5244.036
38	4.5244.038
40	4.5244.040
42	4.5244.042
44	4.5244.044
46	4.5244.046
48	4.5244.048
50	4.5244.050
52	4.5244.052
54	4.5244.054
56	4.5244.056
58	4.5244.058
60	4.5244.060
65	4.5244.065
70	4.5244.070
75	4.5244.075
80	4.5244.080
85	4.5244.085
90	4.5244.090
95	4.5244.095
100	4.5244.100
105	4.5244.105
110	4.5244.110

ChM sp. z o.o.

Lewickie 3b
16-061 Juchnowiec Kościelny
Polska

tel. +48 85 86 86 100

fax +48 85 86 86 101

chm@chm.eu

www.chm.eu



CE 0197