ChM®



SPINE STABILIZATION

- IMPLANTS
- INSTRUMENTS 15.0907.011
- INSTRUMENTS 15.0907.012
- INSTRUMENTS 15.0907.002
- ADDITIONAL INSTRUMENTS
- SURGICAL TECHNIQUE



www.chm.eu

SYMBOLS DESCRIPTION



Caution - pay attention to a special procedure.



Perform the activity under X-Ray control.



Information about the next stages of a procedure.



Proceed to the next stage.



Return to the specified stage and repeat the activity.



Before using the product, carefully read the Instructions for Use. It contains, among others, indications, contraindications, side effects, recommendations and warnings related to the use of the product.



The above description is not a detailed instruction of conduct. The surgeon decides about choosing the operating procedure.

www.chm.eu

 Document No
 ST/96A

 Date of issue
 10.11.2020

 Review date
 P-006-06.11.2025

The manufacturer reserves the right to introduce design changes.

Updated INSTRUCTIONS FOR USE are available at the following website: ifu.chm.eu

1. INTRODUCTION			
2. IMPLANTS	8		
3. INSTRUMENTS	31		
3.1. CONTAINERS ARRANGEMENT	38		
4. SURGICAL TECHNIQUE	40		
4.1. THORACOTOMY	40		
4.2. ANTERIOR THORACOLUMBAR APPROACH	41		
4.3. ANTERIOR RETROPERITONEAL APPROACH	41		
4.4. POSTERIOR APPROACH TO THE THORACOLUMBAR APPROACH	42		
4.5. APPROACH TO POSTERIOR SUPERIOR ILIAC SPINE	42		
4.6. SCREW SELECTION. PREPARATION OF THE SCREW ENTRY POINT	43		
4.7. INSERTION OF SCREWS. POSTERIOR APPROACH	44		
4.8. INTRODUCTION OF POLYAXIAL SCREWS INTO THE PELVIS	62		
4.9. HOOKS INSERTION	64		
4.10. IMPLANT REMOVAL	67		
4.11. CEMENT AUGMENTATION OF SCREWS (OPTIONAL)	68		



1. INTRODUCTION

CHARSPINE2 Thoracolumbar Spinal Stabilization System is the set of universal spinal fixation implants for thoracolumbar and lumbar spine treatment in skeletally mature patients:

- via posterior approach screw fixation from T1 (73) to S2 hook fixation from T1 (73) to L5
- via anterolateral approach screw fixation from T4 (T6) to L4 (L3)

CHARSPINE2 system consists of:

- implants (screws, hooks, connectors, locking elements, staples, and others),
- instruments for implants insertion,
- instructions for use and surgical technique.

INDICATIONS

CHARSPINE2 implants allow for treatment intended for spinal physiological curvature reconstruction by means of appropriate vertebrae reposition.

Indications for use:

- · degenerative disc disease,
- · spondylolistheses,
- fractures and instabilities,
- deformities (e.g. scolioses or kyphoses),
- tumours,
- stenoses,
- · pseudoarthroses,
- · nonunion following the previous procedures.

CONTRAINDICATIONS

Contraindications may be relative and absolute. One should thoroughly consider the selection of an appropriate implant on the basis of comprehensive assessment of patient's health condition. Some conditions such as spinal infection, morbid obesity, mental disease, alcohol or drug addiction, pregnancy, oversensitivity to metals/ foreign bodies, insufficient tissue coverage or open wound in the operative site may reduce the chances of surgery or make the success impossible.



A detailed list of contraindications is presented in instructions for use (IFU) intended for this device.

WARNINGS

Safety and effectiveness of spinal systems based upon pedicle screw fixation have been established only for pathological spinal conditions caused by significant mechanical instability or deformations requiring surgical fixation.

Safety and effectiveness of these systems for any other conditions are unknown.

It is not always possible to achieve positive results in each and every patient. This especially applies to procedures in which other conditions related to patient's state may make it impossible to achieve the positive results.

The final result is greatly influenced by appropriate patient selection and patient's observance of postoperative recommendations. It is proved that smoking hampers the bone union. Patients should be informed about this correlation and warned about the consequences.



A detailed list of warnings, precautions and postoperative recommendations is presented in instructions for use (*IFU*) intended for this device.








Implants from the **ChM CHARSPINE2** spine stabilization system are designed and tested to be used only with the appropriate **ChM** instruments. This surgical technique is intended as a guide only. As with any surgical procedure, the surgeon should be thoroughly trained before the procedure and must take into consideration the particular needs of each patient.

MAIN FEATURES AND BENEFITS

Presented implants and instruments are intended for posterior and anterolateral approach.


The presented range of implants is made of titanium, titanium alloys and cobalt alloy in accordance with ISO 5832 standards.






The screws with a double thread of CHARSPINE2 system are compatible with the instruments 15.0913 for minimally invasive technique (M/S). Should these screws be introduced using percutaneous method (without the use of a guide rod), please refer to the surgical technique No. ST-86, intended for the CHARSPINE2 MIS system.





Polyaxial screws allow for stable angular fixation of the screw head in each direction.

Uniplanar screws combine the features of mediallateral stiffness of monoaxial screws with the mobility of polyaxial screws in the rostral-caudal direction. Polyaxial screws for pelvis allow for the extension of thoracolumbar spine stabilization and fixation in the pelvic bone. The screws offer an increased asymmetrical range of motion in one of the planes, facilitating the rod-to-screw fixation.





2. IMPLANTS

CHARSPINE2 MONOAXIAL SCREW







CHARSPINE2 LOCKING SCREW







3.6160.000	/

	Len	Ti
	25	3.6150.525
	30	3.6150.530
4.0	35	3.6150.535
	40	3.6150.540
	45	3.6150.545
	25	3.6151.525
	30	3.6151.530
4.5	35	3.6151.535
	40	3.6151.540
	45	3.6151.545
	25	3.6152.525
	30	3.6152.530
	35	3.6152.535
5.0	40	3.6152.540
	45	3.6152.545
	50	3.6152.550
	25	3.6153.525
	30	3.6153.530
	35	3.6153.535
5.5	40	3.6153.540
5.5	45	3.6153.545
	50	3.6153.550
	55	
		3.6153.555
	25	3.6154.525
	30	3.6154.530
	35	3.6154.535
6.0	40	3.6154.540
6.0	45	3.6154.545
	50	3.6154.550
	55	3.6154.555
	60	3.6154.560
	65	3.6154.565
	25	3.6155.525
	30	3.6155.530
	35	3.6155.535
	40	3.6155.540
6.5	45	3.6155.545
	50	3.6155.550
	55	3.6155.555
	60	3.6155.560
	65	3.6155.565
	25	3.6156.525
	30	3.6156.530
	35	3.6156.535
	40	3.6156.540
	45	3.6156.545
	50	3.6156.550
7.5	55	3.6156.555
1.5	60	3.6156.560
	65	3.6156.565
	70	3.6156.570
	75	3.6156.575
	80	3 6156 590

	(Len)	(Ti)
	25	3.6157.525
	30	3.6157.530
	35	3.6157.535
	40	3.6157.540
	45	3.6157.545
	50	3.6157.550
	55	3.6157.555
	60	3.6157.560
8.5	65	3.6157.565
	70	3.6157.570
	75	3.6157.575
	80	3.6157.580
	85	3.6157.585
	90	3.6157.590
	95	3.6157.595
	100	3.6157.601
	25	3.6158.525
	30	3.6158.530
	35	3.6158.535
	40	3.6158.540
	45	3.6158.545
	50	3.6158.550
	55	3.6158.555
0.5	60	3.6158.560
9.5	65	3.6158.565
	70	3.6158.570
	75	3.6158.575
	80	3.6158.580
	85	3.6158.585
	90	3.6158.590
	95	3.6158.595
	100	3.6158.601
	25	3.6159.525
	30	3.6159.530
	35	3.6159.535
	40	3.6159.540
	45	3.6159.545
	50	3.6159.550
	55	3.6159.555
10.5	60	3.6159.560
10.5	65	3.6159.565
	70	3.6159.570
	75	3.6159.575
	80	3.6159.580
	85	3.6159.585
	90	3.6159.590
	95	3.6159.595
	100	3.6159.601



3.6156.580

3.6156.585

3.6156.590

80

85

90



CHARSPINE2 MONOAXIAL REDUCTION SCREW







CHARSPINE2 LOCKING SCREW










3.6160.000	/
------------	----------



80

85


90

3.6167.580

3.6167.585

3.6167.590

	(Len)	(Ti)
	25	3.6168.525
	30	3.6168.530
	35	3.6168.535
	40	3.6168.540
	45	3.6168.545
	50	3.6168.550
	55	3.6168.555
	60	3.6168.560
8.5	65	3.6168.565
	70	3.6168.570
	75	3.6168.575
	80	3.6168.580
	85	3.6168.585
	90	3.6168.590
	95	3.6168.595
	100	3.6168.601
	25	3.6169.525
	30	3.6169.530
	35	3.6169.535
	40	3.6169.540
	45	3.6169.545
	50	3.6169.550
	55	3.6169.555
9.5	60	3.6169.560
7.5	65	3.6169.565
	70	3.6169.570
	75	3.6169.575
	80	3.6169.580
	85	3.6169.585
	90	3.6169.590
	95	3.6169.595
	100	3.6169.601
	25	3.6149.525
	30	3.6149.530
	35 40	3.6149.535 3.6149.540
	45	3.6149.545
	50	3.6149.550
	55	3.6149.555
	60	3.6149.560
10.5	65	3.6149.565
	70	3.6149.570
	75	3.6149.575
	80	3.6149.580
	85	3.6149.585
	90	3.6149.590
	95	3.6149.595
	100	3.6149.601





CHARSPINE system 2

CHARSPINE2 POLYAXIAL SCREW





CHARSPINE2 LOCKING SCREW








_	/
· ·	

	(Len)	(Ti)
	25	3.6170.525
	30	3.6170.530
	35	3.6170.535
4.0		3.6170.535
	40	
	45	3.6170.545
	50	3.6170.550
	25	3.6171.525
	30	3.6171.530
4.5	35	3.6171.535
	40	3.6171.540
	45	3.6171.545
	50	3.6171.550
	25	3.6172.525
	30	3.6172.530
	35	3.6172.535
5.0	40	3.6172.540
	45	3.6172.545
	50	3.6172.550
	55	3.6172.555
	25	3.6173.525
	30	3.6173.530
	35	3.6173.535
5.5	40	3.6173.540
5.5	45	3.6173.545
	50	3.6173.550
	55	3.6173.555
	25	3.6174.525
	30	3.6174.530
	35	3.6174.535
	40	3.6174.540
6.0	45	3.6174.545
	50	3.6174.550
	55	3.6174.555
	60	3.6174.560
	65	3.6174.565
	25	3.6175.525
	30	3.6175.530
	35	3.6175.535
	40	3.6175.540
6.5	45	3.6175.545
	50	3.6175.550
	55	3.6175.555
	60	3.6175.560
	65	3.6175.565
	25	3.6176.525
	30	3.6176.530
	35	3.6176.535
	40	3.6176.540
	45	3.6176.545
	50	3.6176.550
	55	3.6176.555
7.5		
	60	3.6176.560
	<u>65</u>	3.6176.565
	70	3.6176.570
		2 (17(577
	75	3.6176.575
	75 80 85	3.6176.575 3.6176.580 3.6176.585

	(Len)	(Ti)
	25	3.6530.525
	30	3.6530.530
	35	3.6530.535
	40	3.6530.540
	45	
		3.6530.545
	50	3.6530.550
	55	3.6530.555
8.5	60	3.6530.560
	65	3.6530.565
	70	3.6530.570
	75	3.6530.575
	80	3.6530.580
	85	3.6530.585
	90	3.6530.590
	95	3.6530.595
	100	3.6530.600
	25	3.6531.525
	30	3.6531.530
	35	3.6531.535
	40	3.6531.540
	45	3.6531.545
	50	3.6531.550
	55	3.6531.555
9.5	60	3.6531.560
	65	3.6531.565
	70	3.6531.570
	75	3.6531.575
	80	3.6531.580
	85	3.6531.585
	90	3.6531.590
	95	3.6531.595
	100	3.6531.600
	25	3.6532.525
	30	3.6532.530
	35	3.6532.535
	40	3.6532.540
	45	3.6532.545
	50	3.6532.550
	55	3.6532.555
10.5	60	3.6532.560
	65	3.6532.565
	70	3.6532.570
	75	3.6532.575
	80	3.6532.580
	85	3.6532.585
	90	3.6532.590
	95	3.6532.595
	100	3.6532.600



3.6176.590

90

CHARSPINE2 POLYAXIAL REDUCTION SCREW















Ø4.0	Ø4.5	Ø5.0	Ø5.5	Ø6.0	Ø6.5	Ø7.5	Ø8.5	Ø9.5	Ø10.5
				Cole	ours				

	Len	Ti
	25	3.6177.525
	30	3.6177.530
4.0	35	3.6177.535
	40	3.6177.540
	45	3.6177.545
	25	3.6178.525
	30	3.6178.530
4.5	35	3.6178.535
	40	3.6178.540
	45	3.6178.545
	20	3.6179.520
	25	3.6179.525
	30	3.6179.530
5.0	35	3.6179.535
	40	3.6179.540
	45	3.6179.545
	50	3.6179.550
	25	3.6180.525
	30	3.6180.530
	35	3.6180.535
5.5	40	3.6180.540
	45	3.6180.545
	50	3.6180.550
	55	3.6180.555
	25	3.6181.525
	30	3.6181.530
	35	3.6181.535
	40	3.6181.540
6.0	45	3.6181.545
	50	3.6181.550
	55	3.6181.555
	60	3.6181.560
	65	3.6181.565
	25	3.6182.525
	30	3.6182.530
	35	3.6182.535
	40	3.6182.540
6.5	45	3.6182.545
	50	3.6182.550
	55	3.6182.555
	60	3.6182.560
	65	3.6182.565
	25	3.6183.525
	30	3.6183.530
	35	3.6183.535
	40	3.6183.540
	45	3.6183.545
	50	3.6183.550
7.5	55	3.6183.555
	60	3.6183.560
	65	3.6183.565
	70	3.6183.570
	75	3.6183.575
	80	3.6183.580
	85	3.6183.585
	90	3.6183.590


	Len	Ti
	25	3.6533.525
	30	3.6533.530
	35	3.6533.535
	40	3.6533.540
	45	3.6533.545
	50	3.6533.550
	55	3.6533.555
	60	3.6533.560
8.5	65	3.6533.565
	70	3.6533.570
	75	3.6533.575
	80	3.6533.580
	85	3.6533.585
	90	3.6533.590
	95	3.6533.595
	100	3.6533.600
	25	3.6534.525
	30	3.6534.530
	35	3.6534.535 3.6534.540
	40	
	45	3.6534.545
	50	3.6534.550
	55	3.6534.555
9.5	60	3.6534.560
	65	3.6534.565
	70	3.6534.570
	75	3.6534.575
	80	3.6534.580
	85	3.6534.585
	90	3.6534.590
	95	3.6534.595
	100	3.6534.600
	25	3.6535.525
	30	3.6535.530
	35	3.6535.535
	40	3.6535.540
	45	3.6535.545
	50	3.6535.550
	55	3.6535.555
10.5	60	3.6535.560
10.5	65	3.6535.565
	70	3.6535.570
	75	3.6535.575
	80	3.6535.580
	85	3.6535.585
	90	3.6535.590
	95	3.6535.595
	100	3.6535.600



CHARSPINE2 POLYAXIAL SCREW FOR PELVIS















3.6160.000	/


	Len	Ti
	25	3.6514.525
	30	3.6514.530
	35	3.6514.535
	40	3.6514.540
	45	3.6514.545
	50	3.6514.550
6.5	55	3.6514.555
6.5	60	3.6514.560
	65	3.6514.565
	70	3.6514.570
	75	3.6514.575
	80	3.6514.580
	85	3.6514.585
	90	3.6514.590
	25	3.6515.525
	30	3.6515.530
	35	3.6515.535
	40	3.6515.540
	45	3.6515.545
	50	3.6515.550
7.5	55	3.6515.555
7.5	60	3.6515.560
	65	3.6515.565
	70	3.6515.570
	75	3.6515.575
	80	3.6515.580
	85	3.6515.585
	90	3.6515.590






CHARSPINE2 UNIPLANAR SCREW










	Len	Ti
	25	3.6184.525
	30	3.6184.530
4.0	35	3.6184.535
7.0	40	3.6184.540
	45	3.6184.545
	25	3.6185.525
	30	3.6185.530
4.5	35	3.6185.535
	40	3.6185.540
	45	3.6185.545
	25	3.6186.525
	30	3.6186.530
	35	3.6186.535
5.0	40	3.6186.540
	45	3.6186.545
	50	3.6186.550
	25	3.6187.525
	30	3.6187.530
	35	3.6187.535
5.5	40	3.6187.540
3.5	45	3.6187.545
	50	3.6187.550
	55	3.6187.555
	25	3.6188.525
	30	3.6188.530
	35	3.6188.535
	40	3.6188.540
6.0	45	3.6188.545
0.0	50	3.6188.550
	55	3.6188.555
	60	3.6188.560
	65	3.6188.565
	25	3.6189.525
	30	3.6189.530
	35	3.6189.535
	40	3.6189.540
6.5	45	3.6189.545
0.5	50	3.6189.550
	55	3.6189.555
	60	3.6189.560
	65	3.6189.565
	25	3.6190.525
	30	3.6190.530
	35	3.6190.535
	40	3.6190.540
	45	3.6190.545
	50	3.6190.550
	55	3.6190.555
7.5	60	3.6190.560
	65	3.6190.565
	70	3.6190.570
	75	3.6190.575
	80	3.6190.580
	85	3.6190.585
	90	3.6190.590
	20	3.0170.370

CHARSPINE2 UNIPLANAR REDUCTION SCREW










(∠')	(Len)	(Ti)
	25	3.6191.525
	30	3.6191.530
4.0	35	3.6191.535
	40	3.6191.540
	45	3.6191.545
	25	3.6192.525
	30	3.6192.530
4.5	35	3.6192.535
	40	3.6192.540
	45	3.6192.545
	25	3.6193.525
	30	3.6193.530
	35	
5.0		3.6193.535
	40	3.6193.540
	45	3.6193.545
	50	3.6193.550
	25	3.6194.525
	30	3.6194.530
	35	3.6194.535
5.5	40	3.6194.540
	45	3.6194.545
	50	3.6194.550
	55	3.6194.555
	25	3.6195.525
	30	3.6195.530
	35	3.6195.535
	40	3.6195.540
6.0	45	3.6195.545
	50	3.6195.550
	55	3.6195.555
	60	3.6195.560
	65	3.6195.565
	25	3.6196.525
		3.6196.530
	30	
	35	3.6196.535
	40	3.6196.540
6.5	45	3.6196.545
	50	3.6196.550
	55	3.6196.555
	60	3.6196.560
	65	3.6196.565
	25	3.6197.525
	30	3.6197.530
	35	3.6197.535
	40	3.6197.540
	45	3.6197.545
	50	3.6197.550
	55	3.6197.555
7.5	60	3.6197.560
	65	3.6197.565
	70	3.6197.570
		3.6197.575
	80	3.6197.580
	85	
		3.6197.585
	90	3.6197.590




ROD 6





Len	Ti	Co
40	3.3246.040	4.3980.040
50	3.3246.050	4.3980.050
60	3.3246.060	4.3980.060
70	3.3246.070	4.3980.070
80	3.3246.080	4.3980.080
90	3.3246.090	4.3980.090
100	3.3246.100	4.3980.100
120	3.3246.120	4.3980.120
140	3.3246.140	4.3980.140
160	3.3246.160	4.3980.160
180	3.3246.180	4.3980.180
200	3.3246.200	4.3980.200
220	3.3246.220	4.3980.220
260	3.3246.260	4.3980.260
300	3.3246.300	4.3980.300
360	3.3246.360	4.3980.360
400	3.3246.400	4.3980.400
460	3.3246.460	4.3980.460
500	3.3246.500	4.3980.500




Len	Ti	Co
40	3.3248.040	4.3249.040
50	3.3248.050	4.3249.050
60	3.3248.060	4.3249.060
70	3.3248.070	4.3249.070
80	3.3248.080	4.3249.080
90	3.3248.090	4.3249.090
100	3.3248.100	4.3249.100
120	3.3248.120	4.3249.120
140	3.3248.140	4.3249.140
160	3.3248.160	4.3249.160
180	3.3248.180	4.3249.180
200	3.3248.200	4.3249.200



CURVED ROD 6








Len	Ті
35	3.6280.035
40	3.6280.040
45	3.6280.045
50	3.6280.050
55	3.6280.055
60	3.6280.060
65	3.6280.065
70	3.6280.070
75	3.6280.075
80	3.6280.080
85	3.6280.085
90	3.6280.090
95	3.6280.095
100	3.6280.100
110	3.6280.110
120	3.6280.120
130	3.6280.130
140	3.6280.140
150	3.6280.150
160	3.6280.160
170	3.6280.170
180	3.6280.180
190	3.6280.190
200	3.6280.200






Len	Ti
35	3.6295.035
40	3.6295.040
45	3.6295.045
50	3.6295.050
55	3.6295.055
60	3.6295.060
65	3.6295.065
70	3.6295.070
75	3.6295.075
80	3.6295.080
85	3.6295.085



ANGLED ROD 6







	Len	Ti
	150	3.3981.155
	200	3.3981.205
50° —	250	3.3981.255
30	300	3.3981.305
	350	3.3981.355
	450	3.3981.455
	150	3.3981.156
	200	3.3981.206
60° —	250	3.3981.256
- 00	300	3.3981.306
	350	3.3981.356
	450	3.3981.456
	150	3.3981.157
	200	3.3981.207
70° —	250	3.3981.257
/0	300	3.3981.307
	350	3.3981.357
	450	3.3981.457






	Len	Co
	150	4.3981.155
	200	4.3981.205
50°	250	4.3981.255
50	300	4.3981.305
	350	4.3981.355
	450	4.3981.455
	150	4.3981.156
	200	4.3981.206
60°	250	4.3981.256
00	300	4.3981.306
	350	4.3981.356
	450	4.3981.456
	150	4.3981.157
70° –	200	4.3981.207
	250	4.3981.257
/ 0	300	4.3981.307
	350	4.3981.357
	450	4.3981.457



CROSSWISE CONNECTOR SOLID (SET)








Len	Ti
14	3.6290.014
15	3.6290.015
16	3.6290.016
17	3.6290.017
18	3.6290.018
19	3.6290.019
20	3.6290.020
21	3.6290.021
22	3.6290.022
23	3.6290.023
24	3.6290.024
25	3.6290.025
26	3.6290.026
27	3.6290.027
28	3.6290.028
29	3.6290.029
30	3.6290.030




CROSSWISE CONNECTOR REGULATED (SET)












	Len	Ti
S	30.5-33	3.3979.030
Μ	33-38.5	3.3979.033
L	38.5-49	3.3979.038
XL	49-71	3.3979.049



	Len	Ti
	22	3.6296.022
	24	3.6296.024
	26	3.6296.026
XS	28	3.6296.028
	30	3.6296.030
	32	3.6296.032
	34	3.6296.034

	Len	Ti	
XS	33-37.5	3.3972.033	

	Len	Ti
S	37.5-40	3.3972.037
Μ	40-45.5	3.3972.040
L	45.5-56.5	3.3972.045
XL	56.5-78	3.3972.056
XXL	78-99	3.3972.078



AXIAL CONNECTOR (SET)







d	DZ	Ti
5	5	3.3970.855
6	5	3.3970.865
6	6	3.3970.866

PARALLEL CONNECTOR (SET)








d		Ti
5	5	3.3970.955
6	5	3.3970.965
6	6	3.3970.966

PARALLEL CONNECTOR (SET)







d	D7	Ti
6	6	3.6294.012




ANGULAR CONNECTOR









	Len	Ti
0°	12	3.6284.012
0	16	3.6284.016
10°	12	3.6285.012
10	16	3.6285.016
30°	12	3.6286.012
30	16	3.6286.016

CHARSPINE2 LOCKING SCREW










3.6160.000

/

LATERAL CONNECTOR









	Len	Ti
	15	3.6281.015
	20	3.6281.020
90°	25	3.6281.025
	30	3.6281.030
	35	3.6281.035
	15	3.6282.015
	20	3.6282.020
75°	25	3.6282.025
	30	3.6282.030
	35	3.6282.035
	15	3.6283.015
	20	3.6283.020
105°	25	3.6283.025
	30	3.6283.030
	35	3.6283.035









CLAMP CROSSWISE CONNECTOR (SET)









LATERAL CONNECTOR







Len	Ti
35	3.6289.035
40	3.6289.040
45	3.6289.045
50	3.6289.050
55	3.6289.055
60	3.6289.060
65	3.6289.065
70	3.6289.070
80	3.6289.080
90	3.6289.090
100	3.6289.100



SPINAL HOOKS



Small		Standard		Large	
Laminar ho	ok				
	3.6266.001		3.6266.002	H	3.6266.003
Laminar ho	ok, narrow blade				
	3.6267.001		3.6267.002	I	3.6267.003
Laminar ho	ok, extended body				
	3.6268.001*		3.6268.002		3.6268.003*
Laminar ho	ok, offset				
	3.6269.001 – right * 3.6269.101 – left *		3.6269.002 – right 3.6269.102 – left		3.6269.003 – right * 3.6269.103 – left *
Laminar ho	ok, angled blade				
	3.6270.001*		3.6270.002		3.6270.003*
Thoracic lar	ninar hook				
			3.6271.002		
Thoracic lar	ninar hook, narrow blade				
			3.6272.002		
Thoracic lar	ninar hook, offset				
			3.6273.002 – Small offset, right; 3.6273.102 – Small offset, left; 3.6274.002 – Large offset, right; 3.6274.102 – Large offset, left;		

^{*} available as additional item



SPINAL HOOKS



Pedicle hook 3.6275.001 3.6275.002 3.6275.003* Transverse process hook	Small		Standard		Large
	Pedicle ho	ok			
Transverse process hook		3.6275.001		3.6275.002	3.6275.003*
	Transverse	process hook			
3.6276.001 – right* 3.6276.002 – right 3.6276.101 – left* 3.6276.102 – left 3.6276.103 – left*		3.6276.001 – right* 3.6276.101 – left*			





The palettes for implants presented below are not offered as sets (they do not include implants).

40.8064.000		Screw	Size	
PALETTE FOR CHARSPINE2 IMPLANTS - SCREWS		diameter	L	No. of sockets
			25	2
			30 35	2
		4.5	40	2 2
			45	2
			50	2
			35 40	2 2
		_	45	2
		5.0 –	50	2
			55	2
			60 30	<u>2</u> 8
			35	8
		5.5	40	8
	Monoaxial screws	-	45 50	8 8
	9 6		55	8
	با ال		35	8
			40	8
		6.0	45 50	8 8
			55	8
			60	8
		_	35	6
			40 45	6
		6.5	50	6
The state of the s			55	6
1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	1		60 35	6
The second secon	→ ∅		40	6
	. Ø .	7.5	45	6
		7.5	50	6
			55 60	6
			50	2
			55	2
		8.5 -	60 70	2 2
			80	2
			90	2
			50	2
			55 60	2 2
		9.5	70	2
			80	2
			90 30	<u>2</u> 4
		-	35	4
	Polyaxial screws	5.5	40	4
	. ,	ے د. <i>د</i>	45	4
			50 55	4
	JI K		35	4
			40	4
	A	6.0	45	4
			50 55	4 4
			60	4
	盡		35	6
	T _	_	40	6
	批	6.5	45 50	6
	1		55	6
			60	6
	₩ <u></u>		35 40	6
	→ Ø ←	_	40 45	6
	. W .	7.5 –	50	6
			55	6
			60	6

40.8119.000 PALETTE SMALL FOR CHARSPINE2 IMPLANTS - SCREWS		Screw diameter	Size L	No. of sockets
			30	6
	Managarial		35	6
	Monoaxial, poly-	5.0	40	6
	axial and unipla- nar screws		45	6
	Tial sciews		50	5
	31 16		30	6
11/42/22/2019			35	6
		5.5	40	6
			45	6
			50	5
			35	6
		_	40	6
	_	6.0	45	6
		0.0	50	6
			55	6
			60	5
			35	6
	Y		40	6
		6.5	45	6
		0.5	50	6
			55	6
			60	5





40.8065.000 Palette for CHARSPINE2 implants - Connectors 1	Implant type	Size	No. of sockets
		L-40	2
		L-50	2
	Rods	L-60	2
		L-70	2
		L-80	2
		L-90 L-100	2
		L-120	4
		L-140	2
	_	L-160	2
		L-180	2
of the little		L-200	2
		L-220	4
		L-260	4
		L-300	2
		L-360	2
		L-460	2
	Rod connectors	L-60	2
	_	L-80	2
		L-100	2
	Locking screws	-	28
	Axial connector	6/6	1
	Parallel connector	6/6	1
	Clamp crosswise connector	-	4



EXCHANGEABLE MODULES – IMPLANT SOCKETS CONFIGURATION

40.8078.000 Exchangeable module 1	Implant type	Size	No. of sockets
	Curved rod _	L-35	1
		L-40	1
		L-45	1
		L-50	1
		L-55	1
	_	L-60	1
		L-65	1
		L-70	1
		L-75	1
		L-80	1
	<u> </u>	L-85	1

40.6795.000 Exchangeable module 5	Implant type	Size	No. of sockets
	Curved rod	L-90	1
		L-95	1
		L-100	1
		L-110	1
		L-120	1
		L-130	1
		L-140	1
		L-150	1
		L-160	1

40.8080.000 Exchangeable module 3	Implant type	Size	No. of sockets
		L-14	1
		L-15	1
	Crosswise connector solid	L-16	1
		L-17	1
		L-18	1
		L-19	1
		L-20	1
	Crosswise connector regulated -	XS	1
	(monoaxial)	S	1
		М	1
		L	1
		XL	1
	Lateral and a star (and anish	XS	1
	Lateral connector (polyaxial)	S	1
		М	1
		L	1
		XL	1



MODULAR PALETTE FOR CHARSPINE2 IMPLANTS - CONNECTORS 2 (STANDARD CONFIGURATION)









40.8081.000 Exchangeable module 4



40.8079.000 Exchangeable module 2

EXCHANGEABLE MODULES – IMPLANT SOCKETS CONFIGURATION

40.8081.000 Exchangeable module 4	Implant type	Size	No. of sockets
	Angular connector	L-12	1
	Arigular confriector	L-16	1
		10° L-12	1
		10° L-16	1
		30° L-12	1
3.6281 0.50 - 3.6282 0.75 - 3.6283 0.75		30° L-16	1
		90° L-15	1
		90° L-20	1
		90° L-25	1
A Company of the Comp	Lateral connector	90° L-30	1
		75° L-15	1
		75° L-20	1
		75° L-25	1
L-15 L-20 L-25 L-30 L-15 L-20 L-25 L-30 L-15 L-20 L-25 L-30		75° L-30	1
		105° L-15	1
		105° L-20	1
		105° L-25	1
		105° L-30	1



 $It\ is\ possible\ to\ change\ the\ configuration\ of\ modules\ included\ into\ palettes\ according\ to\ an\ individual\ order.$

Right

Left

Transverse process hook

3

3



Lnvi				IMPLANTS
40.8077.000 Palette for CHARSPINE2 implants - Hooks		Implant type	Variant	No. of sockets
40.8077.000			Small	6
		Laminar hook	Standard	6
			Large	6
			Small	6
		Laminar hook, narrow blade	Standard	6
			Large	6
		Laminar hook, extended body	-	2
	TE	Landau la Con	Right	2
		Laminar hook, offset -	Left	2
	P	Laminar hook, angled blade	-	2
		Thoracic laminar hook	-	3
	Y	Thoracic laminar hook, narrow blade	-	3
		Thoracic laminar hook,	Right	3
	-	narrow blade	Left	3
		Thoracic laminar hook,	Right	3
	7	large offset	Left	3
		Pedicle hook	Small	2
		redicte floor	Standard	2



3. INSTRUMENTS

Instruments for CHARSPINE2 spine stabilizer - basic [15.0907.011]					
Instruments for CHARSPINE2 module 1 [15.0907.111]	Name	Catalogue no.	Pcs		
	Container lid 9x4	14.0907.113	1		
	Container 9x4H	14.0907.116	1		
	Screw persuader Screw persuader is used to press the rod down to the bottom of the transpedicular screw cut-out.	40.8096.100	1		
	Thoracic pedicular trocar Thoracic pedicular trocar is used to prepare openings in the pedicle of the vertebral arch in the thoracic section of the spine.	40.8070.000	1		
	Universal pedicular trocar Universal pedicular trocar is used to prepare openings in the pedicle of the vertebral arch in the lumbar section of the spine.	40.8071.000	1		
	Straight pedicular trocar Straight pedicular trocar is used to prepare openings in the pedicle of the vertebral arch in the lumbar section of the spine.	40.8072.000	1		
30 3/0	Holding forceps Holding forceps are used to conduct the rod derotation procedure.	40.6202.000	1		
	Wrench for monoaxial screws Wrench for monoaxial screws is used for insertion and mounting of CHARSPINE2 monoaxial transpedicular screws. It is intended for use with T-type or oval head ratchet handle.	40.6734.000	1		
	Wrench for polyaxial screws Wrench for polyaxial screws is used for insertion and mounting of CHARSPINE2 polyaxial transpedicular screws. It is intended for use with T-type or oval head ratchet handle.	40.6737.000	2		
	Adjustable rod bender Adjustable rod bender is used to bend the rod to desired shape.	40.8074.000	1		
	Oval head ratchet handle Oval head ratchet handle is used with wrenches for screws and cortical tap.	40.8086.000	2		



Instruments for CHARSPINE2 spine stabilizer - basic [15.0907.011]			
Instruments for CHARSPINE2 module 2 [15.0907.112]	Name	Catalogue no.	Pcs
	Container 9x4H	14.0907.112	1
	Counter wrench Counter wrench is used to ensure rotational stability of the implants system during final tightening of the locking screws.	40.8095.000	1
	Screwdriver tip T30 Screwdriver tip T30 is intended to be used with T-type torque handle 12Nm [40.8087.000]. It is used to finally lock the transpedicular screws, hooks and lateral connectors.	40.8084.000	1
	Screwdriver T30 Screwdriver T30 is used for application and initial locking of the locking screws.	40.8111.000	1
	Pedicle probe straight Pedicle probe is used to verify the continuity of the vertebral arch pedicle.	40.6698.000	1
~**************************************	Pedicle probe curved Pedicle probe is used to verify the continuity of the vertebral arch pedicle.	40.6699.000	1
	Compression forceps-jaws W-26 (set) Exchangeable compression jaws are used with compression forceps.	40.5768.026	1
	Compression forceps-jaws W-46 (set) Exchangeable compression jaws are used with compression forceps.	40.5768.046	1
	Distraction forceps-jaws Exchangeable distraction jaws are used with distraction forceps.	40.5769.000	1
	Parallel distraction forceps Distraction forceps are used with exchangeable jaws and are intended for procedure of vertebrae distraction.	40.8093.000	1
	Parallel compression forceps Compression forceps are used with exchangeable jaws and are intended for procedure of vertebrae compression.	40.8094.000	1
	Trocar Trocar is used to puncture the cortical layer of the vertebral arch pedicle, as a point of insertion of transpedicular screw.	40.8073.000	1



Instruments for CHARSPINE2 spine stabilizer - basic [15.0907.011]							
Instruments for CHARSPINE2 module 4 [15.0907.211]	Name	Catalogue no.	Pcs				
Lanca I	Tray 4x4 1/2H	14.0907.211	1				
	Pliers for rod Pliers for rod are used to grab and insert the spinal rod.	40.8109.000	1				
	Torque limiting ratchet handle T 12Nm Torque limiting ratchet handle T 12Nm is intended to be used with screwdriver tip T30 [40.8084.000] and is used for final tightening of the locking screws into the transpedicular screws, hooks and lateral connectors.	40.8087.000	1				

Instruments for CHARSPINE2 spine stabilizer - extension [15.0907.012]			
Instruments for CHARSPINE2 module 3 [15.0907.114]	Name	Catalogue no.	Pcs
	Container lid 9x4	14.0907.115	1
	Container 9x4H	14.0907.114	1
	Cortical tap 4.0	40.8075.540	1
	Cortical tap 4.5	40.8075.545	1
	Cortical tap 5.0	40.8075.550	1
***************************************	Cortical tap 5.5	40.8075.555	1
***************************************	Cortical tap 6.0	40.8075.560	1
	Cortical tap 6.5	40.8075.565	1
<	Cortical tap 7.5	40.8075.575	1
	Cortical tap 8.5	40.8075.585	1
	Cortical tap 9.5 Cortical taps are intended for usage with ratchet handles [40.8085.000] or [40.8086.000] and may be used for tapping the vertebral arch pedicle prior to screw insertion.	40.8075.595	1
FOR LOS CONTACTOR SOUR MODIFICACION CONTRACTOR (CONTRACTOR CONTRACTOR CONTRAC	Screwdriver tip T15 Screwdriver tip T15 is intended to be used only with T-type torque handle 3.5Nm [40.8088.000]. It is used to screw the locking screws into crosswise, axial and parallel connectors.	40.8110.000	1
	Reduction screw device Reduction screw device is used to break off the extended arms of reduction screws.	40.8108.000	1
This is a second of the second	Rod bender 6.0 left Bender is used for bending the rod in situ.	40.8091.000	1



Instruments for CHARSPINE2 spine stabilizer - extension [15.0907.012]				
Instruments for CHARSPINE2 module 3 [15.0907.114]	Name	Catalogue no.	Pcs	
CAN	Rod bender 6.0 right Bender is used for bending the rod in situ.	40.8092.000	1	
	Fork persuader Fork persuader is used to press the rod down to the bottom of the transpedicular screw cut-out.	40.8100.000	1	
	Eye wrench Eye wrench is used to prevent the rod from changing its position while bending, and to conduct the rod derotation.	40.8069.000	1	
	Rod impactor Rod impactor is used to impact and press the rod down to the transpedicular screw cut-out.	40.8068.000	1	
Instruments for CHARSPINE2 module 5 [15.0907.212]	Name	Catalogue no.	Pcs	
The state of the s	Tray 4x4 1/2H	14.0907.212	1	
	Holder for crosswise connectors Holder for crosswise connectors is used for insertion and mounting of clamp and regulated crosswise connectors.	40.8067.000	1	
	Holder for rod connectors Holder for rod connectors is used for insertion and mounting of rod connectors (which are mounted together with clamp connectors).	40.8076.000	1	
	T-type torque handle 3.5Nm T-type torque handle 3.5Nm is intended to be used only with screwdriver tip T15 [40.8110.000] and is used for final tightening of the locking screws into the crosswise, axial and parallel connectors.	40.8088.000	1	
EWW TEMPLATE A SERVICE SERVIC	Template Template is used for selecting the size of crosswise and rod connectors.	40.5248.000	1	



Instruments for CHARSPINE2 spine stabilizer - extension 1 [15.0907.002] (instruments for Spinal hooks) Instruments for CHARSPINE2 module 6 [15.0907.203] Name Catalogue no. Pcs Tray 5x4 1/2H Hook holder Hook holder Hook holder is used to insert spinal hooks. Lateral hook holder Lateral hook holder is used to insert spinal hooks. 40.8102.000 1

Instruments for CHARSPINE2 module 7 [15.0907.204]	Name	Catalogue no.	Pcs
	Tray 5x4 1/2H	14.0907.204	1
	Impactor for hooks Impactor for hooks is used for final impaction of spinal hook into the selected space.	40.8103.000	1
	Narrow raspatory for laminar hooks Narrow raspatory for laminar hooks is used to prepare space for a laminar hook.	40.8104.000	1
	Raspatory for laminar hooks Raspatory for laminar hooks is used to prepare space for a laminar hook.	40.8105.000	1
	Wide raspatory for laminar hooks Wide raspatory for laminar hooks is used to prepare space for a laminar hook.	40.8106.000	1
	Raspatory for pedicle hooks Raspatory for pedicle hooks is used to prepare space for a pedicular hook.	40.8107.000	1



Extended instruments [15.0907.002] is additional equipment.

In order to include the instruments to the ordered **CHARSPINE2** instruments, please contact your local representative or **ChM** Sales Department.



Instruments mentioned below are not included in the standard instrument set.
In order to include them to the ordered **CHARSPINE2** instruments, please contact your local representative or **ChM** Sales Department.

Supple	ementary instruments	
	Name	Catalogue no.
	Screw persuader Screw persuader may be used as an alternative to the standard persuader 40.8096.100. The instrument can be handled using only one hand.	40.8083.100
. F. F. T. 134 - F. T. 1594 - F. F. F. T. 1594 - F. F. F. T. 1594 - F. F. F. F. T. 1594 - F.	Rod trial 6/300 Rod trials are used for initial rough assessment of the size and shape of the rod and to facilitate the selection of the proper size of the spinal rod, in the spinal stabilization procedures using transpedicular screws.	40.5246.300
	Pedicle probe Pedicle probe is used to verify the continuity of the vertebral arch pedicle. The probe is equipped with two tips of different stiffness.	40.6696.000
	Tubular rod bender (2 pcs.) Tubular rod bender is used for bending the rod to the acute angle, especially when stabilization: the sacral bone – iliac ala is performed.	40.6178.000
	Distraction forceps Szczypce dystrakcyjne są używane do przeprowadzenia procedury dystrakcji kręgów.	40.6176.000
	Compression forceps Compression forceps are used for vertebrae compression.	40.6694.000
	Screwdriver T30 short Screwdriver T30 short is used as an alternative to the standard screwdriver T30 [40.8111] in situations where the conditions of surgery or the surgeon's preferences require the use of a shorter instrument.	40.6151.100
	Wrench for monoaxial screws short Short wrench for monoaxial screws is used for insertion and mounting of CHARSPINE2 monoaxial transpedicular screws. It is intended for use with T-type or oval head ratchet handle in space-restricted conditions.	40.8920.000
	Wrench for polyaxial screws short Short wrench for polyaxial screws is used for insertion and mounting of CHARSPINE2 polyaxial transpedicular screws. It is intended for use with T-type or oval head ratchet handle in space-restricted conditions.	40.8921.000



Supple	ementary instruments	
	Name	Catalogue no.
	Thoracic pedicular trocar Thoracic pedicular trocar [40.6243] is used as an alternative to the standard thoracic trocar [40.8070] in situations where the surgeon's preferences require the use of an instrument with an oval handle.	40.6243.000
	Universal pedicular trocar Universal pedicular trocar [40.6244] is used as an alternative to the standard universal trocar [40.8071] in situations where the surgeon's preferences require the use of an instrument with an oval handle.	40.6244.000
	Straight pedicular trocar Straight pedicular trocar [40.6245] is used as an alternative to the standard straight trocar [40.8072] in situations where the surgeon's preferences require the use of an instrument with an oval handle.	40.6245.000
	Hand hold rod cutter Hand hold rod cutter is used for easy cutting of rods with diameters of 6mm, 5mm and 3.5mm.	40.5288.000
	Container lid 9x4	14.0907.106
	Container 9x4H Use the container [14.0907.105] with silicone mat [12.0725.000] for storage of the additional instruments (provided on the customer's request) of the CHARSPINE2 system. This container, depending on the needs, can be included into another container of the CHARSPINE2 system or provided with a separate container lid [14.0907.106].	14.0907.105
	Perforated aluminum lid	12.0750.200
	Container with solid bottom Container used for sterilization of instruments [15.0907.012]	12.0750.101
	Container with solid bottom Container used for sterilization of instruments [15.0907.011]	12.0750.103

Additional instruments, not included in the set, for bone cement delivery through an open access

Name Catalogue no. Pcs



In order to include them to the ordered CHARSPINE2 instruments, please contact your local representative or ChM Sales Department.


	Cannula for bone cement		
	Single use cannula for bone cement is intended to be used with CHARSPINE2 MIS screws. The universal Luer thread allows the cannula to be connected with a mixing and delivery system for bone cement.	40.8591.000	1
	Alignment trial for cannula		
	Alignment trial for cannula, inserted through the cannula for bone cement attached to the CHARSPINE2 MIS screw, is intended to confirm the coaxiality of the holes in the cannula and screw.	40.8592.000	1
N	Pusher	40.8596.000	1
	Pusher is used to remove residual bone cement from the cannula.	40.6596.000	,
	Counter wrench		
	Counter wrench is used to immobilize the head of the polyaxial screw when installing/ removing the cannula for bone cement.	40.6749.000	1



3.1. CONTAINERS ARRANGEMENT


Instruments - basic [15.0907.011]

No.	Name	Catalogue No.	
1	Container lid 9x4	14.0907.113	1
2	Container 9x4H	14.0907.116	1
3	Tray 4x4 1/2H	14.0907.211	1
4	Container 9x4H	14.0907.112	1



Instruments extended [15.0907.012]

No.	Name	Catalogue No.	Pcs
1	Container lid 9x4	14.0907.115	1
2	Tray 4x4 1/2H	14.0907.212	1
3	Container 9x4H	14.0907.114	1






Instruments extended [15.0907.002]

Option I				
No.	Name	Catalogue No.	Pcs	
1	Instruments - basic	15.0907.011	1	
2	Tray 5x4 1/2H	14.0907.203	1	
3	Instruments - extended	15.0907.012	1	
4	Tray 5x4 1/2H	14.0907.204	1	



Option II (In case extended Instruments 15.0907.012 has not been purchased)				
No.	Name	Catalogue No.	Pcs	
1	Perforated aluminum lid	12.0750.200	1	
2	Tray 5x4 1/2H	14.0907.203	1	
3	Tray 5x4 1/2H	14.0907.204	1	
4	Container with solid bottom	12.0750.100	1	





4. SURGICAL TECHNIQUE

Anterior approach to thoracolumbar spine

Surgical procedures on the thoracolumbar spine by means of anterior approach are generally performed with a patient in a lateral position, with the assistance of a general or vascular surgeon.

4.1. THORACOTOMY

Thoracotomy is a standard approach for the treatment of thoracic spine disorders such as deformity, tumor or infection. In case of deformity treatment, the approach is always located on the side of the curve apex, e.g. a right-sided thoracotomy is chosen for a right-sided curve. In general, a left-sided thoracotomy is preferred, especially in the lower thoracic area, due to right-sided location of the liver which limits the operative field. However, when the upper part of the thoracic spine is concerned, some surgeons favour right-sided approach (*in cases when the spinal pathology does not dictate the side of thoracotomy*) to avoid subclavian and carotid arteries in the left superior mediastinum.


Indications


Spinal pathologies (deformities, degenerations, fractures, tumours, infections) that are located between T4 and T10 are indications for the thoracotomy.

Patient positioning

In case of right-sided thoracotomy, the patient is lying on the left side on a soft, rubber mattress. The arms are positioned at elevation of 90° and with the elbows flexed. The legs are kept straight, with the right leg resting on the left leg. The symphysis and the sacrum are supported by pads to maintain the specified position. Prior to skin incision, the side of thoracotomy and the level involved are to be confirmed. It is essential to center the incision right over the pathology place and to select the intercostal space correctly. To confirm the selected spine level, it is recommended to count the ribs and to compare the result with the radiograph.

The skin incision shall be extended from the lateral border of the paraspinous muscle up to the sternocostal joint.







4.2. ANTERIOR THORACOLUMBAR APPROACH

The anterior approach to thoracolumbar section may be used if there is a need of simultaneous exposure of vertebral bodies of lower thoracic and upper lumbar parts of the spine. Technically, this approach is more difficult than thoracotomy because of the diaphragm exposed and the increased risk of simultaneous exposure of the thoracic cavity and the peritoneal space. If the spine pathology does not determine the side of the approach, the access from the left side is preferred due to right-sided location of the liver.

Indications

The anterior thoracolumbar approach is recommended for spine pathologies mentioned as an indication for thoracotomy and situated between T9 and L5.


Patient positioning

The patient is placed in the right lateral decubitus position, with supports placed beneath the thorax and shoulders. The table may be slightly bent above the level of pelvis to increase the distance between pelvis and thorax.

During the operation special care should be taken to not harm the branches of the phrenic nerve, which are extending peripherally from the center towards anterolateral and posterior direction. It is recommended to make the incision around the periphery of the diaphragm to minimize the interference with its function when making the thoracoabdominal approach to the spine.

Special care should also be taken when entering the abdominal cavity.

To gain the best access to the space between T12 and L1, it is usually recommended to resect the tenth rib which allows exposure between T10 and L2.



4.3. ANTERIOR RETROPERITONEAL APPROACH

The anterior retroperitoneal approach to the lumbar vertebral bodies is a modification of the anterolateral approach commonly used by general surgeons during the sympathectomy. It allows for superior, multilevel access to the lumbar spine.

Indications


The anterior retroperitoneal approach is recommended for spine pathologies (*deformities, degenerations, fractures, tumours, infections*) situated between L2 and L5.


Patient positioning

The patient is placed in the decubitus position, in most cases on the right side. Most often, the approach is made from the left side to prevent damage to the liver and the inferior vena cava. To better expose the space between the twelfth rib and the iliac crest, the table planes may be flexed. Lower limbs are bent slightly in hips to release the tension of the psoas muscle.

The incision is to be oblique, above the twelfth rib, from the lateral border of the quadratus lumborum muscle to the lateral border of the rectus abdominis muscle, in order to allow access to the first and second lumbar vertebrae.

When the lower vertebrae (from L3 to L5) are exposed, the incision is to be made a few fingers below and parallel to the costal margin.







4.4. POSTERIOR APPROACH TO THE THORACOLUMBAR APPROACH

The posterior approach to the thoracolumbar spine can be made through standard midline longitudinal incision with lateral retraction of the erector spinae in the direction of the transverse processes tips. This approach allows for access to the spinous processes, vertebral arches and joints at all levels.

The target spine level should be determined using the X-Ray control, so that the spine is unveiled only at the required segment.


Indications


The anterior approach to the thoracolumbar spine is recommended for spine pathologies (*deformities, spine canal stenoses, fractures, degenerations, tumours, infections, instabilities, herniations*) situated between T1 and L5.

Patient positioning

Patient is placed prone on rubber-foam supports. To avoid excessive pressure and pressure sores, a headrest with support for mouth, nose and eyes should be used. It is vital to avoid any pressure on the abdomen. It is crucial while decompressing the spine, as pressure on abdomen may cause vein congestion and thus excessive intraoperative bleeding.

Positioning the patient on a bending surgical table with supports with flexion of hip and knee joints allows for reduction of lumbar lordosis and easier access to posterior spine elements and spine canal structures, especially at the lumbosacral junction.







4.5. APPROACH TO POSTERIOR SUPERIOR ILIAC SPINE

Indications

This approach is recommended when the following occur: a significant lumbopelvic instability (caused by damage at S1 level resulting from trauma, tumour or infection) or long thoracolumbosacral instrumentation of scoliosis, causing a high risk of instability of the lumbosacral connection.

Patient positioning

Patient is positioned in the same manner as presented in section 4.4.

Screw implantation in pelvis requires access to the posterior superior iliac spine. First, the lumbosacral spine is exposed. The posterior superior iliac spine may be exposed with a separate, longitudinal skin incision, bilateral resection of the myofascial flaps and retraction in lateral and cephalad direction.

The entry point is located at the lower part of the posterior superior iliac spine. It is recommended to use osteotome (*or rongeur*) to remove a fragment of the iliac crest around the screw head or to sink the screw head in the bone to avoid any screw prominence, especially when slim patients are concerned.



4.6. SCREW SELECTION. PREPARATION OF THE SCREW ENTRY POINT

During transpedicular stabilization it is of vital importance to select appropriate screw diameter for specific vertebrae and to carefully choose the site and α angle of insertion. Depending on the location level, the pedicles of vertebrae arches are varied in terms of shapes and geometry (e.g. the cross section of the vertebrae arch pedicles in the thoracic spine indicates an irregular, kidney-like shape with the medially-directed convexity).


Taking into consideration the above-mentioned, the initial selection of screw diameter and length has to be performed within the preoperative procedures, individually for each vertebrae on the basis of CT and X-Ray images (in AP and lateral projections).

The internal dimension of the arch of a vertebrae pedicle (*W*) is of vital importance when choosing the external diameter of the transpedicular screw. It is crucial to remember that the pedicle dimensions obtained on the basis of imaging in AP projection are not real dimensions and should be treated as approximate values only. In general, the outer diameter of the screw is 2 mm smaller than the internal dimension of the vertebrae pedicle arch.

Screw insertion point is located at the intersection of a line that divides the transverse processes in half and the line along the lateral aspect of the superior articular process.



The surgeon shall decide about the size of screws on the basis of CT and X-Ray imaging and intraoperative identification (*probing the pedicle*).



There are two alternative trajectories for insertion of screws through the thoracic vertebrae pedicles:

A – anatomical approach


B – straight approach (direct)

The insertion point is located at the intersection of a line in sagittal projection about 1 mm in medial direction from the lateral edge of the lamina and of a line along the transverse processes about 1 mm below the surface of the superior articular process.



If anatomical approach is used, only polyaxial screws are to be used.

If straight (*direct*) approach is used, both normal and polyaxial screws may be used.





4.7. INSERTION OF SCREWS. POSTERIOR APPROACH

4.7.1. PREPARATION OF VERTEBRAL ARCH PEDICLES




The point of screw insertion is prepared with a trocar [40.8073.000] which is used to puncture the cortical layer of the vertebral arch pedicle.

When it is necessary a bone rongeur is used to remove the upper part of the vertebral articular process at the screw insertion point, therefore the cancellous bone right under the cortical layer and the access to the vertebral arch pedicle are exposed.

Pedicle diameter and the angle should be determined prior to the operation by means of imaging studies. It allows for later determination of depth and angle of the prepared canal and the screw diameter.





An opening for screw is prepared with the use of a pedicular trocar (which is available as: universal [40.8071.000]/[40.6244.000], straight [40.8072.000]/ [40.6245.000] and thoracic [40.8070.000]/[40.6243.000]).



Trocars [40.6243.000], [40.6244.000], [40.6245.000] are additional equipment and are not part of the instruments for the CHARSPINE 2 system.

The instrument is inserted by means of delicate rotary-oscillatory movement.

The tip should be inserted carefully, led along the interior walls of the vertebral cortical bone with the smallest resistance possible, so the vertebra walls remain undamaged.



Trocar tip has marked depth indicators in five-millimeter increments to help determine the correct length of the transpedicular screw.

Identical procedure should be used while preparing the opening in the second pedicle.










Prior to screw insertion, it is recommended to check the continuity of all walls of the vertebral arch pedicle with the help of the pedicle probe [40.6698], [40.6699] or [40.6696.000].




Pedicle probe [40.6696.000] constitutes additional equipment and is not included in the instruments for the CHARSPINE 2 system.










CHARSPINE2 screws are self-tapping, therefore there is no need to tap the vertebral arch pedicle. However, if tapping is clinically required, it is possible to use taps [40.8075.540]÷[40.8075.595] mounted on oval head ratchet handle [40.8086.000].




Before threading, make sure that the tap is compatible with the screw that will be implanted.





4.7.2. INSERTION OF SCREWS





Wrenches for monoaxial and polyaxial screws are intended to be mounted on:

• oval head ratchet handle [40.8086.000].



Monoaxial and polyaxial wrenches have a mechanism that prevents any spontaneous loosening of the tip-screw connection during the transpedicular screws insertion.



Wrenches [40.8920.000] and [40.8921.000] constitute additional equipment and are not included in the instruments for the CHARSPINE 2 system.

Wrench for monoaxial screws [40.6734.000] / [40.8920.000]

Wrench for polyaxial screws [40.6737.000] / [40.8921.000]



Square end of the wrench is mounted in the quick coupling end of the handle **[40.8086.000]**.



Then an appropriate length and diameter of the transpedicular screw is selected.

The tip is inserted all the way into the screw head:

• in the case of monoaxial screws, a tip of wrench for monoxial screws is to be used.




• in the case of polyaxial and uniplanar screws, a tip of wrench for polyaxial screws is to be used.






By turning the knob clockwise, continue tightening the threaded outer sleeve until the screws is set still in relation to the wrench.



The red ring on the sleeve of the wrench should be covered by the screw. Should the ring be not covered by the screw, loosen the wrench by turning the knob counterclockwise, re-position the wrench tip in the screw socket and then retighten the wrench.



Then pull the slider to secure the screw from spontaneous loosening.





The screw mounted on a wrench is inserted into an opening prepared beforehand.




Screw insertion should be controlled in two planes with the help of X-Ray control.



Precise screw positioning should be performed by tightening the screws, not by loosening them.



Moving back the screw may result in loss of connection stability and may necessitate the use of a larger diameter screw.






The core of the transpedicular screw is strengthened in the vicinity of its head.

To reduce the potential risk of screw breakage, it is necessarry to screw it all the way in so the whole thread is in the bone.


To uninstall the wrench from the screw,

- 1. press the slider and then use counterclockwise
- 2. rotation to unscrew the threaded sleeve of the wrench.



CAUTION:

In case of difficulties when pressing the slider:



- 2. Hold the ratchet handle to fix the transpedicular screw, then rotate slightly the knob under the slider in a clockwise direction (this will relieve load on the directional coupling).
- 3. Press the slider and unscrew the threaded sleeve of the wrench.



4.7.3. ROD SHAPING





Benders 40.6178.000 constitute additional equipment and are not included in the instruments for the CHARSPINE 2 system.


Having inserted the screws, select a rod of appropriate length.




In order to determine the approximate length and the desired shape of the rod, rod trial [40.5246.300] can be used.

To achieve the desired spine curvature (e.g. lordosis or kyphosis), the rod should be appropriately shaped. Shaping is performed with the help of adjustable rod bender [40.8074.000] or tubular benders [40.6178.000].


To secure the rod against movement during shaping, eye wrench [40.8069.000] should be used.














Appropriately shaped rods are inserted into cut-outs of transpedicular screws with the help of pliers for rod [40.8109.000].





To correct the shape of the rod in situ, use the rod benders - right [40.8092.000] and left [40.8091.000].





If necessary, cut the rod to the desired length with the use of hand held rod cutter [40.5288].



Hand hold rod cutter is a non-standard instrument and is not included into the **CHARSPINE2** Instruments.



4.7.4. ROD FIXATION

The rod is locked by inserting the locking screw [3.6160.000] into the transpedicular screw head.




The locking screw may be mounted on the screwdriver tip only from the upper side of the screw (the locking screw design eliminates any errors related to the mounting).



The upper surface of the screw is colour-marked to allow for easier identification.






The locking screw is mounted on the tip of the screwdriver T30 [40.8111.000] or screwdriver [40.6151.000], then it is inserted into the cut-out on the screw head and screwed in slightly in a clockwise direction, simultaneously gently pressing the rod to the screw cutout bottom.



Screwdriver [40.6151.000] constitutes additional equipment and is not included in the instruments for the CHARSPINE 2 system.






In case of difficulties when pressing the rod to the screw cut-out bottom, it is possible to use rod impactor [40.8068.000].












To press the rod to the screw cutout bottom, it is possible to use also fork persuader [40.8100.000].

The fork persuader [40.8100.000] is of special design that allows its use even when the adjacent screws are very close to each other.








- the rod must completely adhere to the cut-out bottom in the screw head,
- the upper part of the locking screw (in blue) should flush with the upper part of the screw head.




It is necessary to avoid the following:

- the rod is not placed horizontally in the screw head,
- the rod is high and does not adhere to the bottom of the screw head cut-out,
- the screw is embedded in the place of rod bending (on the convexity or concavity of the arch).




When the osteosynthesis covers more than three vertebrae, there is a risk that the rod may not fit to all screws cut-outs. In such case reduction screws are of help (with prolonged, breakable head arms). Use of a few reduction screws facilitates correct positioning and locking of the rod.






When the long arms spread while inserting the locking screw, first put the counter wrench [40.8095.000] on the screw head and press the rod, then screw in the locking screw until the rod is completely pressed to the screw head bottom.





Screws arms are broken off at the end of the surgery with the use of reduction screw device [40.8108.000].






The abovementioned procedure may cause excessive reduction of screws (and vertebrae); to avoid this it is necessary to correct the rod bending in situ.



The use of polyaxial screws allows for adjusting the screw head position in relation to the rod, especially when the screws are not parallel. In this manner the connection achieves greater stability and there is no need for complicated rod bending. It is possible to tilt the screw head in any direction.

The use of uniplanar screws gives the possibility to tilt the head of the screw in the rostral-caudal direction, while ensuring medial-lateral stiffness.







It is possible to use the holding forceps **[40.6202.000]** during the procedure of rod derotation.

In such case, having established the desired rod position, the rod should be locked to maintain its position.

This allows for the next stage to be performed - the reposition of vertebrae.







At this stage it is possible to perform:

• the vertebrae distraction with the use of parallel distraction forceps [40.8093.000] or distraction forceps [40.6176.000],


• the vertebrae compression with the use of parallel compression forceps [40.8094.000] or compression forceps [40.6694.000].



A proper shaping of rods is a crucial stage of surgery allowing for good vertebrae reposition.



The compression forceps [40.6694.000] and distraction forceps [40.6176.000] constitute additional equipment and are not included in the instruments for the CHARSPINE 2 system.






If more force is necessary to impact the rod into the bottom of the transpedicular screw cut-out, it is possible to use screw persuader [40.8096.100].



Prior to use, extend the clamping arms of the persuader as far as possible by turning the knob counterclockwise until a distinct resistance is felt.




Incompletely extended clamping arms of the persuader may cause the instrument to be improperly connected with the screw and may damage the device.



Then, by turning clockwise, the rod may be smoothly pushed down to the bottom of the screw.






40.8111.000

Then turn the persuader knob counterclockwise to loosen the device and extend the clamping arms completely.

To secure the rod, insert the locking screw through the cannulated opening of the screw persuader (*the screw is mounted on the screwdriver tip T30* [40.8111.000]) and initially screw it in.




The device is dismounted from the transpedicular screw by skewing the device in the rostral-caudal direction.



Incompletely extended clamping arms of the persuader may damage the device while dismounting.





Having established the required vertebrae position, finally screw in the locking screws with the help of T-type torque handle 12Nm [40.8087.000] connected with screwdriver tip T30 [40.8084.000].

When the required torque of 12Nm is reached, the dynamometric mechanism signals it with an audible snap.

To eliminate the movement of rod-screws system while screwing in the locking elements use the counter wrench [40.8095.000].




To maintain high level of safety and correctness of operation of the torque wrench, it is necessary to keep the calibration date presented on the stopper of the instrument handle.

The instrument calibration is performed by the manufacturer - **ChM** sp. z o.o.




If it is needed to lengthen the fixation in lateral direction (in relation to the main axis of stabilization), it is possible to use a lateral connector. The connector is put on the main rod, then it is locked in a desired position (after mounting the appropriate transpedicular screw).



When it is necessary to use two sections of rods (e.g. in case of scolioses), it is possible to use the following connectors:


- crosswise connector,
- parallel connector.





CHARSPINE2 allows for increase of rotational stability of the system by connecting two sides of rods with the help of a crosswise connector.

CHARSPINE2 offers three types of crosswise connectors:

- 1. regulated crosswise connectors.
- **2.** clamp crosswise connectors with rod connectors.
- **3.** solid crosswise connectors (*intended for dual-rod stabilization with an anterior approach*).











40.5248.000

To facilitate the selection of appropriate length of the crosswise connector, use the template [40.5248.000]




4.8. INTRODUCTION OF POLYAXIAL SCREWS INTO THE PELVIS

Polyaxial screws for pelvis extend the lumbosacral stabilization and provide the fixation in the pelvic bone. The screw offers an increased, asymmetrical range of motion in one of the planes, facilitating screw-to-rod fixation.

The screw insertion hole is prepared using a universal [40.8071.000] or straight pedicular trocar [40.8072.000].



Pedicular trocar may not provide an opening corresponding to the full length of the screw for pelvis. This should be confirmed intraoperatively by X-Ray imaging.





The optimum insertion trajectory of the screw for pelvis starts above the incisura ischiadica major and ends in the quadrilateral area of the pelvic bone above the linea glutea inferior. This trajectory ensures optimum fixation of the screw in the bone.





Polyaxial screw of a determined size is screwed into the prepared pelvic opening using the wrench for polyaxial screws [40.6737.000], attached to the oval handle [40.8086].


It is recommended to use osteotome (*or rongeur*) to remove a fragment of the iliac crest around the screw head or to sink the screw head in the bone to avoid any screw prominence, especially when slim patients are concerned.






Most often the standard straight rod, e.g. [3.3246.xxx] is used for stabilization. The rod should be profiled so that the screw inserted into the S1 vertebra can be connected with the screw inserted into the pelvis.

It is also possible to use factory profiled rod, e.g. [3.3981.xxx] should any difficulties with rod shaping appear.



A lateral connector, e.g. [3.6283.xxx] may also be used to connect the rod with the screw inserted into the pelvis.

Locking of the polyaxial screw for pelvis is performed in the same way as for standard polyaxial screws.





4.9. HOOKS INSERTION

4.9.1. PEDICLE HOOKS INSERTION

Pedicle hook is used for fixation of thoracic vertebrae and it is recommended for spine section between T10 and T1 vertebrae. Pedicle hook is always inserted in rostral direction, so its split blade leans on the vertebral arch pedicle.

Place for pedicle hook insertion is prepared by a limited facetectomy. On the desired level, there are two incisions made on the surface of the inferior transverse process, enabling the access to the cartilage of the superior transverse process of the preceding vertebra.




Hook insertion point may be prepared with the use of raspatory for pedicle hooks [40.8107.000] by careful insertion in slightly lateral direction (*in relation to the medial line*) until the pedicle is identified. It is vital not to penetrate medially the spine canal.





When the place of insertion is ready, the hook is inserted with the help of hook holder [40.8101.000], and then it is carefully impacted to the desired position with the use of impactor for hooks [40.8103.000].








4.9.2. INSERTION OF LAMINAR HOOKS

Laminar hooks are used in thoracolumbar spine. They may be inserted in rostral or caudal direction, depending on the spine section instrumented. A wide range of laminar hooks is available. The selection of appropriate hooks depends on the anatomy of the insertion point:

- in case of a hook inserted in caudal direction (*supralaminar manner*), it is recommended to use the thoracic hook with narrow blade to avoid excessive penetration of the spine canal with the hook blade.
- offset laminar hooks are recommended in situations when standard hooks do not ensure collinearity of inserted implants.
- extended laminar hooks are used in situation when a specific height (in relation to other implants) has to be maintained.

Laminar hook blade is inserted in space above the dura mater. To allow for appropriate passage of hook in spinal canal, a ligamentum flavum is removed and a limited laminectomy is performed.





During the preparation of space for hooks the following raspatories for laminar hooks may be used:

- narrow [40.8104.000],
- standard [40.8105.000],
- wide [40.8106.000].

Blade widths correspond to widths of laminar hooks available.





The selected hook is mounted in jaws of a hook holder [40.8101.000] and then implanted in a prepared site in vertebral pedicle.






Alternatively, the selected hook may be mounted in jaws of lateral hook holder [40.8102.000], e.g. in situation when anatomical structure hinders the use of a holder [40.8101.000].

4.9.3. INSERTION OF TRANSVERSE PROCESS HOOKS

Transverse process hooks are often used in thoracic spine region because of fairly large size of transverse process. Hooks may be inserted in rostral or caudal position. When the hook is inserted on a transverse process in caudal direction, it may be used for insertion (*in one line*) with inferiorly inserted pedicle hook to achieve an appropriate fixation and to ensure better stability.

Therefore, between the transverse process and the rib in front of the process, an upper and front surface of the transverse process is prepared with the use of raspatory for laminar hooks [40.8105.000].






Mount the selected hook on a hook holder [40.8101.000] and then insert it on a prepared transverse process.



4.10. IMPLANT REMOVAL

To perform revision, the following steps should be taken (observing the order provided).

- 1. Use the screwdriver T30 [40.8111.000] to loosen and remove the locking screws.
- **2.** Use the pliers for rod [40.8109.000] to remove the rods.
- **3.** Remove the transpedicular screws. For this purpose, a wrench for monoaxial screws **[40.6734.000]** or wrench for polyaxial screws **[40.6737.000]** should be used.





4.11. CEMENT AUGMENTATION OF SCREWS (OPTIONAL)



If bone cement is to be injected into the vertebral body, use fenestrated screws offered in a separate ChM implants system - CHARSPINE2 MIS (catalog pages for fenestrated screws are provided in a separate surgical technique - ST-86).


The CHARSPINE2 MIS fenestrated screws are installed to the wrench 40.6737.000 in the same way as the CHARSPINE2 system screws are.

For the instruments below to be supplied, please contact your local representative or ChM Sales Department.

4.11.1. USE OF CANNULA FOR BONE CEMENT For a bone cement to be injected into the vertebral body, install the cannula for bone cement [40.8591.000] to the implanted screw. To do this, hold the flattened tip of the cannula and insert the other end, into the screw head and then tightened it up (clockwise rotation). While installing the cannula, use the counter wrench [40.6749.000] to immobilize the screw head. 40.8591.000 40.6749.000

To confirm the coaxiality of the holes in the cannula and the screw, use the alignment trial for cannula [40.8592.000] (not included in the standard instrument set). The trial should be inserted into the cannula until stop.











Should the screw been inserted so deep that the screw head was immobilized by the neighbouring bone, it may be necessary to slightly back out the screw to allow the head to move and the cannula to be properly inserted.



Incorrect connection of the cannula with the screw can cause cement leakage at the connection point. Only use **ChM** cannula for bone cement supply.

Before cement injection, all cannulas for bone cement should be attached to the screws.



4.11.2. CEMENT PREPARATION AND INJECTION

For cement preparation, refer to the Instructions for Use for bone cement and cement mixing/delivery device system. The cannula for bone cement is equipped with a standardized Luer Lock thread, enabling a tight connection with the cement mixing/delivery device.



The volume of cement in the cannula [40.8594.000] is 1.2 ml.

Mix the cement as instructed and suck it into the cement mixing/delivery device. Before injecting, wait until the cement reaches the right viscosity.



Fluoroscopy should be used throughout the cement injection procedure to control the cement flow.

When using cement to augment multiple screws and levels, attention must be paid not to exceed the working time of the cement prior to completion of cement delivery through the screws. When the cement working time is close to complete, a new cement package should be used.



The cannula for bone cement as well as the mixing/delivery device are designed for use with one package of bone cement only. If a second package of cement is needed, use a new cannula and mixing/delivery device.

Do not attempt to force the injection of cement if excessive resistance is felt. Always determine the cause of the resistance and take appropriate action. If the cement is seen outside of the vertebral body or in the circulatory system during the procedure, immediately stop injecting the cement.



Use pusher [40.8596.000] to remove the residual cement from the cannula $(1.2\ ml)$. Insert the pusher down the cannula. The pusher handle needs to rest against the cannula.



CAUTION:

Wipe clean the pusher after each use, thoroughly.





4.11.3. REMOVAL OF THE CANNULA FOR BONE CEMENT


After cementing, unscrew the cannula from the screw. Hold the counter wrench [40.6749.000] with the other hand to counteract the removal process.

Back out the cement feeder slightly to ensure that after cementing and before removing the cannula from the screw, the cement flow has been stopped.



It is critical that no torsion movement should be applied to the screw after injecting the cement in order to avoid breaking the cement bridges between the screw and bone.

The cannula for bone cement and cement mixing/delivery device are disposable equipment and must be discarded after cementing.



ChM sp. z o.o.

Lewickie 3b 16-061 Juchnowiec Kościelny Poland tel. +48 85 86 86 100 fax +48 85 86 86 101 chm@chm.eu www.chm.eu



C € ₀₁₉₇