

ChM produces and distributes advanced medical solutions in 3 main divisions:



Comprehensive portfolio of products for traumatology, including systems designed for fractures fixation and deformities correction of extremities and pelvis.



Wide range of advanced solutions for cervical and thoraco-lumbar stabilization of spine, including pedicle screw systems for open and MIS procedures, various interbody devices and fixation plates.



Instruments and implants for cranio-maxillofacial surgeries, dedicated for fracture fixations, reconstructions, distractions and orthognathic surgeries.



ACJ PLATE

3.7242.xxx



To bring medical solutions

ACJ PLATE



ACJ PLATE

Design aligned with anatomy

- optimized anatomical profile
- fit to a wide range of anatomies
- facilitated anatomical reduction

Multiple plate options

- 3 hook heights are offered for perfect fit to the acromioclavicular joint anatomy

Plate clamp

- facilitated minimal surgical access implantation
- reliable stabilization at the trauma site

Stabilizing spikes

- secure anchoring
- stable implant position

Chamfered plate borders

- minimized soft-tissue irritation
- improved stress distribution

Beveled tip

- easy percutaneous insertion

Locking hole

- additional implant stabilization

Set screw


- reliable stabilization at the trauma site
- low head profile to reduce soft tissue irritation



Design aligned with anatomy



Multiple plate options



Stabilizing spikes



Set screw