


# CHM<sup>®</sup>

Instrument set for removing CHARFIX/CHARFIX2 nails  
40.5850.500



[www.chm.eu](http://www.chm.eu)

DO NOT STERILIZE THIS CHART

						Instruments					
						1	2	3	4		
Nail name	Cat. no. of locking elements (1- steel, 3- titanium)	Diameter or size	Locking element socket			Screwdriver tip or wrench			Connector	Impactor- extractor	Mallet
			S	T	Spec.	S	T	Spec.			
<b>ChFN2</b> Trochanteric nail	1.5804.xxx/3.5804.xxx	10.5	-	-	Spec. .	-	-	40.5053.000	40.4731.000	40.5507.000	40.3667.000
	1.5961.xxx /3.5961.xxx	10.5									
	1.5805.xxx /3.5805.xxx	5.0									
	1.5159.xxx /3.5159.xxx	5.0									
	1.5161.003 /3.5161.003	M8	-	25	-	-	-	40.5985.000			
	1.5962.000/3.5962.000	M4									
	1.5161.6xx/3.5161.6xx	M12x1.75									
	1.5808.000/3.5808.000	M6									
Locking mechanism	M12	7	-	-	40.5964.000	-	-				
<b>ChFN</b> Trochanteric nail	1.1949.xxx/3.1949.xxx	11	-	-	Spec. .	-	-	40.5053.000	40.4731.000	40.5507.000	40.3667.000
	1.1938.xxx/3.1938.xxx	11									
	1.1948.xxx/3.1948.xxx	6.5	4	-	-	40.5983.000	-	-			
	1.1935.xxx/3.1935.xxx	6.5									
	1.1657.xxx/3.1657.xxx	5.0	3.5	-	-	40.5982.000	-	-			
	1.1654.xxx/3.1654.xxx	4.5									
	1.2104.6xx/3.2104.6xx	M12x1.75	4	-	-	40.5983.000	-	-			
	1.2106.008/3.2106.008	M8x1.25									
1.2104.003/3.2104.003	M8x1.25	3.5	-	-	40.5982.000	-	-				
<b>CHARFIX2</b> Anatomical femoral nail	1.5168.xxx/3.5168.xxx	6.5	-	30	-	-	40.5986.000	-	40.5512.000	40.5507.000	40.3667.000
	1.5160.xxx/3.5160.xxx	5.5									
	1.5159.xxx/3.5159.xxx	5.0	-	25	-	-	40.5985.000	-			
	1.5161.7xx/3.5161.7xx	M10x1.5	-	30	-	-	40.5986.000	-			
	1.5162.000/3.5162.000	M10x1.5									
<b>CHARFIX2</b> Condylar retrograde femoral nail	1.5151.xxx/3.5151.xxx	6.5	-	25	-	-	40.5985.000	-	40.5857.000	40.5507.000	40.3667.000
	1.5160.xxx/3.5160.xxx	5.5									
	1.5159.xxx/3.5159.xxx	5.0									
	1.5172.000/ 3.5172.000	-	10	-	-	40.5965.000	-	-			
	1.5158.xxx/3.5158.xxx	6.5	-	25	-	-	40.5985.000	-			
	1.5161.006/3.5161.006	M8 spec.									
<b>CHARFIX</b> Femoral nail	1.1652.xxx/3.1652.xxx	6.5	5	-	-	40.5984.000	-	-	40.5071.000	40.5507.000	40.3667.000
	1.2109.xxx/3.2109.xxx	6.5									
	1.1651.xxx/3.1651.xxx	6.5									
	1.1653.xxx/3.1653.xxx	4.5	3.5	-	-	40.5982.000	-	-			
	1.1654.xxx/3.1654.xxx	4.5									
	1.2104.002/3.2104.002	M10x1									
	1.2104.3xx/3.2104.3xx	M10x1	5	-	-	40.5984.000	-	-			
	1.2106.007/3.2106.007	M10x1									

						Instruments					
						1			2	3	4
Nail name	Cat. no. of locking elements (1- steel, 3- titanium)	Diameter or size	Locking element socket			Screwdriver tip or wrench			Connector	Impactor- extractor	Mallet
			S	T	Spec.	S	T	Spec.			
<b>CHARFIX2</b> Retrograde tibial nail	1.5160.xxx/3.5160.xxx	5.5									
	1.5159.xxx/3.5159.xxx	5.0		25		-	40.5985.000	-	40.5857.000		
	1.5161.006/3.5161.006	M8 spec.									
	1.5162.006/3.5162.006	M7x1									
<b>CHARFIX2</b> Tibial nail	1.5160.xxx/3.5160.xxx	5.5									
	1.5159.xxx/3.5159.xxx	5.0									
	1.5170.xxx/3.5170.xxx	4.5	-	25	-	-	40.5985.000	-	40.4710.000		
	1.5169.xxx/3.5169.xxx	4.0									
	1.5161.xxx/3.5161.xxx	M8x1.25									
	1.5162.002/3.5162.002	M8x1.25									
<b>CHARFIX</b> Tibial nail and retrograde tibial nail	1.1657.xxx/3.1657.xxx	5.0									
	1.1654.xxx/3.1654.xxx	4.5									
	1.2104.1xx/3.2104.1xx	M8x1.25	3.5	-	-	40.5982.000	-	-	40.4710.000		
	1.2106.002/3.2106.002	M8x1.25									
	1.2104.003/3.2104.003	M8x1.25									
<b>CHARFIX2</b> Humeral nail	1.5170.xxx/3.5170.xxx	4.5									
	1.5169.xxx/3.5169.xxx	4.0									
	1.5173.xxx/3.5173.xxx	3.0	-	25	-	-	40.5985.000	-	40.5858.000		
	1.5161.2xx/3.5161.2xx	M6x1									
	1.5162.004/3.5162.004	M6									
<b>CHARFIX</b> Humeral nail	1.1657.xxx/3.1657.xxx	5.0									
	1.1654.xxx/3.1654.xxx	4.5									
	1.1653.xxx/3.1653.xxx	4.5									
	1.1655.xxx/3.1655.xxx	3.5	3.5	-	-	40.5982.000	-	-	40.5029.000		
	1.2104.5xx/3.2104.5xx	M7									
	1.2104.2xx/3.2104.2xx	M7									
	1.2106.004/3.2106.004	M7x1									
<b>CHARFIX</b> Forearm and clavicle nail	1.1203.xxx/3.1203.xxx	2.7									
	1.1022.xxx/3.1022.xxx	1.5/2.7									
	1.2104.004/3.2104.004	M4x0.7	2.5	-	-	40.5981.000	-	-	40.5051.000		
	1.2106.006/3.2106.006	M4									

Document No ST-516A  
Date of issue 12.01.2021  
Review date P-001-12.01.2021

*The manufacturer reserves the right to introduce design changes.*

**ChM sp. z o.o.**  
Lewickie 3b  
16-061 Juchnowiec Kościelny  
Polska  
tel. +48 85 86 86 100  
fax +48 85 86 86 101  
chm@chm.eu  
www.chm.eu

