

ChM produces and distributes advanced medical solutions in 3 main divisions:



Comprehensive portfolio of products for traumatology, including systems designed for fractures fixation and deformities correction of extremities and pelvis.



Wide range of advanced solutions for cervical and thoraco-lumbar stabilization of spine, including pedicle screw systems for open and MIS procedures, various interbody devices and fixation plates.



Instruments and implants for cranio-maxillofacial surgeries, dedicated for fracture fixations, reconstructions, distractions and orthognathic surgeries.

## ANGULAR SCREWDRIVER

40.6436.000



To bring medical solutions

[www.chm.eu](http://www.chm.eu)

Document No B/N-01  
Review date P-001-19.08.2020



## ANGULAR SCREWDRIVER

ChM Micro Plates  
**ChMP** system

# ANGULAR SCREWDRIVER

## Increased durability

- perfectly matched bearings minimize wear
- optimized shaft for torque transmission
- improved gear design

## Reduced head profile

- only 8mm head height
- operating flexibility
- minimally invasive

## Fastening clamp

- fast exchange of drill bits and screwdriver tips
- ease of use

## Holding fork

- facilitated implants handling and positioning
- removable

## Rotation handle

- perfect sense of tightening torque
- ergonomic shape

## Connection to a power source


- fast and easy connection
- compatible with ISO 3964

## Minimized number of elements

- easy to clean
- easy to assemble and disassemble

## Comprehensive instrumentation

- full range of complementary devices
- color-coded complementary devices
- transparent module for drill bits and screwdriver tips
- storage tray



Fastening clamp/Reduced head profile

Holding fork

Connection to a power source

Rotation handle