

ChM produces and distributes advanced medical solutions in 3 main divisions:



Comprehensive portfolio of products for traumatology, including systems designed for fractures fixation and deformities correction of extremities and pelvis.



Wide range of advanced solutions for cervical and thoraco-lumbar stabilization of spine, including pedicle screw systems for open and MIS procedures, various interbody devices and fixation plates.



Instruments and implants for cranio-maxillofacial surgeries, dedicated for fracture fixations, reconstructions, distractions and orthognathic surgeries.

ANGULAR SCREWDRIVER

40.6436.000



To bring medical solutions

www.chm.eu

Document No B/N-01
Review date P-001-05.08.2020



ANGULAR SCREWDRIVER

ChM Micro Plates
ChMP system

ANGULAR SCREWDRIVER

Increased durability

- perfectly matched bearings minimize wear
- optimized shaft for torque transmission
- improved gear design

Reduced head profile

- only 8mm head height
- operating flexibility
- minimally invasive

Fastening clamp

- fast exchange of drill bits and screwdriver tips
- ease of use

Holding fork

- facilitate implants handling and positioning
- removable if necessary

Rotation handle

- perfect sense of tightening torque
- ergonomic shape

Connection to a power source


- fast and easy connection
- compatible to standard ISO 3964

Minimized number of elements

- easy to clean
- easy to assemble and disassemble

Comprehensive instrumentation

- full range of complementary devices
- complementary devices with color coding
- transparent module for drill bits and screwdriver blades
- storage tray



Fastening clamp/Reduced head profile

Holding fork

Connection to a power source

Rotation handle