

ChM produces and distributes advanced medical solutions in 3 main divisions:



Comprehensive portfolio of products for traumatology, including systems designed for fractures fixation and deformities correction of extremities and pelvis.



Wide range of advanced solutions for cervical and thoraco-lumbar stabilization of spine, including pedicle screw systems for open and MIS procedures, various interbody devices and fixation plates.



Instruments and implants for cranio-maxillofacial surgeries, dedicated for fracture fixations, reconstructions, distractions and orthognathic surgeries.



ENDOSTEAL PLATES

3.7060.xxx; 3.7061.xxx



To bring medical solutions

www.chm.eu

Document No B/P-08
Review date P-001-04.08.2020

ENDOSTEAL PLATES



ENDOSTEAL PLATES

Design aligned with anatomy

- optimized anatomical profile
- fit to a wide range of anatomies
- facilitated anatomical reduction

Multiple plate options

- system consists of 2 types of plates:
 - curved version
 - straight version
- different variants of length

Curved fixation profile

- stable fixation in the intramedullary canal
- increased implant stiffness

Longitudinal rib

- improved fixation in the intramedullary canal
- angular rotation stability

Sharp tip

- facilitated insertion into the intramedullary canal

2 anatomically directed screws

- stable fixation of in the metatarsal head
- avoided collisions with articular surface
- enhanced plate and fixation stability

Special instruments

- specially designed for surgical technique
- dedicated targeter, raspator and other devices

Multiple screw options

- locking and VA locking screw gives multiple configuration for individual cases

Variable-Angle screws

- high strength cobalt alloy material
- compatible with all locking holes
- 30° angulation cone
- VA screw re-lock possibility



Design aligned with anatomy

Plate-stabilizing rib

Sharp tip

Variable-Angle screws