

ChM produces and distributes advanced medical solutions in 3 main divisions:



Comprehensive portfolio of products for traumatology, including systems designed for fractures fixation and deformities correction of extremities and pelvis.




Wide range of advanced solutions for cervical and thoraco-lumbar stabilization of spine, including pedicle screw systems for open and MIS procedures, various interbody devices and fixation plates.



Instruments and implants for cranio-maxillofacial surgeries, dedicated for fracture fixations, reconstructions, distractions and orthognathic surgeries.

DISTAL LATERAL FIBULA PLATES

3.4158.xxx



To bring medical solutions



DISTAL LATERAL FIBULA PLATES

DISTAL LATERAL FIBULA PLATES

Design aligned with anatomy

- optimized anatomical profile
- fit to a wide range of anatomies
- facilitated anatomical reduction

Chamfered plate borders

- minimized soft-tissue irritation

Beveled tip



- easy percutaneous insertion

Compression hole

- oblong hole for plate positioning
- compression in two directions possible

K-wire holes

- provisional plate positioning



Design aligned with anatomy

Compression hole

K-wires hole